Manual plastic drencher

FOLLOW THE INSTRUCTIONS
Please follow the care and maintenance instructions carefully. This is the only way to ensure the proper functioning of your drencher syringes over time. Keep this manual in a safe place for future reference.

INTENDED USE
Manual drencher syringes are typically used to administer liquid treatments, such as pesticides, vaccines, or nutritional supplements, to cattle, horses, goats, sheep, and other large animals. The drencher syringe is a widely used tool in veterinary treatments that allows liquid solutions to be administered directly into the animal's oral cavity.

CHARACTERISTICS
The drencher syringe has an ergonomic design to ensure the animal's comfort and can be disassembled for easy cleaning.
Drencher syringe components:



1. Plastic dosage adjustment ring
2. Stainless steel spout
3. Plastic graduated cylinder
4. Plastic handle
5. Silicone gasket

HOW TO USE
1. Set the desired dosage using the adjustment ring on the handle; the dosage scale is visible on the front of the cylinder.
2. Check the recommended dosage of the liquid to be administered and ensure that the drencher syringe is adjusted to the prescribed dose..
3. Press the handle and immerse the tip of the nozzle in the liquid to be administered.
4. Release the handle and the integrated spring will draw in the liquid to be administered up to the preset dosage.
5. Make sure that the amount sucked in does not exceed the set dosage.
6. At this point the syringe is ready for use.

AFTER USE
Clean the drencher syringe immediately after use with common detergents (a soft cloth or cotton wool can be used) and rinse thoroughly with warm water. This is essential to avoid residue buildup or contamination with previously administered medications and to ensure proper functioning of all components the next time you use them.

ANIMAL SAFETY TIPS
Hold the animal still; you may ask for the support of a second person to keep the animal calm and still.
Make sure the animal is comfortable throughout the process to avoid stress or injury.

PRECAUTIONS
1. Read the label before unpacking the product.
2. Do not use the drencher syringe if it is damaged.
3. The use of the drencher syringe is reserved for adequately trained and qualified professionals.
4. Before performing any operation, wear protective gloves and wash your hands with antibacterial and antimicrobial soap.
5. Avoid using acidic cleaners, solvents, mineral oils, WD-40, glycerin, or alcohol. Petroleum-based oils can also damage the seal and reduce the product's lifespan.
6. Do not sterilize the syringe. If necessary, you can sterilize the tip.
7. Do not leave the drencher syringe wet for a prolonged period of time.
8. Do not autoclave this product.
9. If leaks occur, check gaskets and connections for signs of wear or cracks.
10. Avoid exposure to direct sunlight.

LUBRICATION
1. Lubricate the seal if the drencher syringe tends to jam.
2. The drencher syringe can be disassembled and each part can be cleaned and lubricated easily.
3. Unscrew the cylinder from the handle and remove it from the plunger.
4. Apply silicone oil to the gasket.
5. Place the cylinder back onto the plunger and screw it into the handle.

LATEX FREE
The product and accessories are latex-free.

NO REPAIR
This product contains no end-user serviceable components.

DISPOSAL
This product must be disposed of in accordance with waste disposal regulations and waste management policies.

GIMA WARRANTY TERMS
The Gima 12-month standard B2B warranty applies.

Note:
This is not a medical device, but a device intended for veterinary use only.

EXPLANATION OF SYMBOLS:

Symbol	Description	Symbol	Description
	Follow the instructions for use		Non sterile
	Intended use		See the instructions for use
	Characteristics		Keep away from sunlight
	Instructions for use		Keep your tools dry
	Safety Tips	EC REP	Authorized Representative in the EU
	Warnings / Precautions	LOT	Lot number
	Lubricant	REF	Catalog Number / Reference
	Latex free		Manufacturer name and address
	No repairs		Manufacturing date
	Disposal	CC	Manufacturing Country

Made in Pakistan

Manual plastic drencher
300 ml



Manual plastic drencher
150 ml



Gima 80170 - 80172
M80170-EN-REV.0.03.26

www.gimaitaly.com

Imported by:
Gima S.p.A.
Via Marconi, 1 - 20060 Gessate (Mi) - Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

9/306 Rang Pura Road, Sialkot 51310, Pakistan
Tel: +92 52 3300220, TEL: +92 333 8636240
Email: info@rizbro.com, Website: www.rizbro.com

REF RA-01-205-15
RA-01-205-30

EC REP CS Active
Usmas iela,19 dz 24 Riga
LV-1083 Lettonia