Gima S.p.A.
Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

PROFESSIONAL MEDICAL PRODUCTS

COMBAT APPLICATION TOURNIQUET

User's manual



1 Band loop

2 Windlass rod

3 Rod locking clip م

4 Self adhering band

5 Rod securing strap -

Notice: When applied in accordance with direction,this tourniquet is a safe and effective device for controling life-threatening extremity bleeding. If you cannot be sure or cannot take the additional time to examine where the bleeding is coming from based on the situation, the tourniquet can be effectively applied over clothing as hight on the arm or leg as possible. The tourniquet must not be applied over solid objects within the clothing. As soon as the situation permits,the injured limb should be evaluated and the tourniquet re-positioned 2-3 inch above the injury directly to the skin.

بلص مسج يأ بلع قبصاعلا عضو متي الأ بجي .قاسلا وأ عارذلا قوف تاصوب 3و 2 نيب حوارتي دعب بلع قبصاعلا .قرشابم دلجلا بلعو قباصملا



1. Apply tourniquet proximal to the bleeding site and place the wounded limb through the band loop.

2. Pull the band very tight and fasten it back on itself.





3. Adhere the band around the limb.



4. Twist the rod until bright red bleeding stops.



5. Place the rod inside the clip and lock.



6. Adhere the band over the rod,for small extremitives,continue to adhere the band around the extremity and over the windlass rod.

7. Secure the rod and band with strap,adhere it to the opposite hook on the windlass clip. Record time of application.



	Conservare in luogo fresco ed asciutto. Keep in a cool, dry place. Á conserver dans un endroit frais et sec. Conservar en un lugar fresco y seco. Armazenar em local fresco e seco. An einem kühlen und trockenen Ort lagern. Διατηρείται σε δροσερό και στεγνό περιβάλλον يحفظ في مكان بارد وجاف	CE	Prodotto conforme alla Direttiva Europea. Product complies with European Directive. Ce produit est conforme à la directive européenne . Producto conforme a la Directiva Europea. Produto em conformidade com a Diretiva Europeia. Produkt gemäß EU-Richtlinie. Προϊόν σύμφωνο με την Ευρωπαϊκή Οδηγία. المنتج متوافق مع التوجيهات الأوروبية
	Leggere e seguire attentamente le istruzioni per l'uso. Read instructions carefully. Lire et suivre attentivement la notice. Leer y seguir atentamente las instrucciones de uso. Ler e seguir atentamente as instruções de uso. Gebrauchsanweisung sorgfältig lesen und befolgen. Διαβάστε και ακολουθήστε προσεχτικά τις οδηγίες χρήσης اقرأ بدقة وحرص تعليمات الاستخدام		Fabbricante. Manufacturer. Fabricant. Fabricante. Fabricante. Hersteller. Παραγωγός. الشركة المصنعة
	Conservare al riparo dalla luce solare. Keep away from sunlight. Á conserver à l'abri de la lumière du soleil. Conservar al amparo de la luz solar. Guardar ao abrigo da luz solar. Vor Sonneneinstrahlung geschützt lagern. Κρατήστε το μακριά από ηλιακή ακτινοβολία يحفظ بعيدًا عن ضوء الشمس		Leggere attentamente le istruzioni per l'uso. Please read instructions carefully. Lire attentivement la notice. Leer atentamente las instrucciones de uso. Ler atentamente as instruções de uso. Gebrauchsanweisung sorgfältig lesen. Διαβάστε προσεχτικά τις οδηγίες χρήσης اقرأ بدقة وحرص تعليمات الاستخدام
REF	Codice prodotto. Product code. Code produit. Código producto. Código produto. Erzeugniscode. Κωδικός προϊόντος كود المنتج	LOT	Numero di lotto. Lot number. Numéro de lot. Número de lote. Número de lote. Chargennummer. Αριθμός παρτίδας رقم الدفعة (انظر العبوة / المغلف)

REF **25526 - 25527**

Gima S.p.A. - Via Marconi, 1
20060 - Gessate (MI) - Italy
Made in China

M25526 - 25527 - M - Rev.0.01.19