

ISO SERVICE TROLLEYS



Size: 65x47xh 100 cm
Weight: 12.5 kg

Pivoting castors
Ø 100 mm
(2 with brake)

27935
MADE IN ITALY



Size: 59x99.3xh 88 cm

Adaptable to any model of ISO 600x400 baskets and trays thanks to the special rail sides, which allow an easy sliding and pull out.

27934

4 wheels Ø 125 mm
(3 with brake)

Lockable with double key



27936

Equipped with key lock and doors kit made of transparent polycarbonate.

Side plastic runners for "ISO" trays or baskets

(baskets and trays not included)

Handle for transport

27936



27924

27938

27938

27938

27939

Antistatic castors
Ø 125 mm (2 with brake)

• **27935 ISO SERVICE TROLLEY - empty**
Service trolley for "ISO" trays or baskets (not included) made of epoxy powder painted steel. Supplied with 5 height adjustable runners for ISO trays/baskets.



Both side basket pulling out with stopper to avoid accidental fall.

• **27934 AURION MINI 600 - ISO trolley - empty**
Basic body made of Baydur Polyurethane - Bayer AG: it can be washed and disinfected. It does not react with diluted acids, alkali and solvents. Self-extinguishing according to class UL V0 - Latex Free. For safety, handles are moulded into the body on both sides, trolley top has a 4-sided moulded containing rim to avoid any spill or fall of materials. No sharp edges. All around bumper. Lightweight and easy to move.


• **27936 ISO TRANSPORT TROLLEY - empty**
Transport trolley for ISO trays and baskets (not included). Made of stainless steel AISI 304 sheet, scotch brite finishing. Anti-dust gasket. Internal side with plastic runners for ISO trays or baskets. Supplied with handle for transport. External size: 52x68xh 159 cm
Internal size: 40x60xh 129 cm

ISO FORMAT UNIVERSAL MODULAR TRAY SYSTEM

27996

ISO MODULAR TRAYS GREY ABS		
Closed	Open	SIZE mm
27924	-	600x400x50
27925	27938	600x400x100
27928	27939	600x400x200
TRANSPARENT PC		
Closed	Open	SIZE mm
27957	-	600x400x50
27958	27965	600x400x100
-	27966	600x400x200

Intelligent solutions for storing all kind of drugs, bandages, accessories and small components in a space saving and clear manner.



27925

27965

27938

GIMA code	DIVIDERS
27994	400x50 mm
27995	600x50 mm
27996	400x100 mm
27997	600x100 mm
27998	400x200 mm
27999	600x200 mm
28001	Drawer stopper

NEO EMERGENCY TROLLEY

• **45722 NEO EMERGENCY TROLLEY**
Heavy duty construction emergency trolley size 75x47.5x92 cm with five drawers. Four castors (Ø 10 cm) assembled on the vertical column line for better stability. One piece counter top with a fully guarded handle for comfort and protection.



Not included

Panel height 8 cm
Inner: 43x33.5x6.8 cm
Panel height 12 cm
Inner: 43x33.5x11 cm

- STANDARD ACCESSORIES**
- Defibrillator plate
 - Sliding side shelf
 - Sharp container
 - Oxygen bottle holder
 - Power outlet
 - Utility container
 - IV pole
 - CPR board
 - 2 waste bins
 - Dividers
 - Central lock



45722

Panel height 24 cm
Inner: 43x33.5x22 cm

NEO MEDICINE TROLLEY

• **45725 NEO MEDICINE TROLLEY - double side**
Two sided medicine trolley 75x48x92 cm with steel and aluminium structure inside, ABS and antibacterial material outside. 4 nylon castors Ø 12.5 cm, one with brake.



Central locking

Aluminium alloy handle

10 drawers (5 per side) with dividers to create up to 20 compartments per drawer



1 large drawer with 12 compartments

45725

- STANDARD ACCESSORIES**
- Central lock
 - 2 waste bin
 - Utility container
 - 4 covered castors
 - Dividers