

HOT AIR STERILIZERS: GIMETTE 1.5, 7, 21 AND 28 LITRES



GIMETTE 1,5

35601



GIMETTE 7

35607

- 35601 GIMETTE 1.5 - hot air sterilizer
- 35607 GIMETTE 7 - hot air sterilizer

Automatic Gimette are compact and high quality hot air sterilizers with the same reliability and efficiency as larger models. Outside structure is made of painted plate with anti-scratch epoxy powder. Stainless steel lid and sterilization chamber with approx. capacity of 1.5 l and 7 l. Removable, perforated for air circulation, tray with pliers. All components certified by most important European brands, assembling in compliance with EN-61010. A second safety thermostat triggers in case temperature should exceed 220°C. Red LED indicates any anomaly in thermoregulation and consequent starting of safety thermostat.

- 35631 RESISTOR 250 W for 35601 - spare
- 35632 RESISTOR 400 W for 35607 - spare
- 35612 GASKET for 35607 - spare
- 35613 GASKET for 35601 - spare

TECHNICAL SPECIFICATIONS

GIMETTE 1.5

Inside dimensions: 245x120x50 mm = vol. 1.5 dm³
 Outside dimensions: 310x180x180 mm Weight: 4.3 kg
 Power consumption: 250 Watt Fuse: 1 A

GIMETTE 7

Inside dimensions: 410x180x90 mm = vol. 7 dm³
 Outside dimensions: 490x310x180 mm Weight: 8 kg
 Power consumption: 400 Watt Fuses: 2 x 2.5 A
 Thermometer: high precision with capillary bulb

Same specifications for both models:

Operating Voltage: 230 V - 50/60 Hz
 Thermostat: automatic, adjustable 0-200 °C
 Timer: manual from 0 to 120 minutes
 2° Thermostat: bi-metallic certified
 Green LED: timer, orange LED: heating, red LED: alarm
 User/Service manual in IT, GB, FR



GIMETTE 21

35621



GIMETTE 28

35628

- 35621 GIMETTE 21 - hot air sterilizer
- 35628 GIMETTE 28 - hot air sterilizer

Fully automatic hot air sterilizer, with adjustable temperature from 0°C to 200°C, kills off even particularly resistant micro-organisms. Stainless steel structure is powder painted, except for the door, to assure maximum heat and shock resistance over time. Four perforated, anodized aluminium trays make hot air circulation easier. Capacity approximately 20 l and 28 l. A second safety thermostat triggers in case temperature should exceed 230°C.



Red LED indicates any anomaly in thermoregulation and consequent starting of safety thermostat. User/Service manual in IT, GB, FR.

TECHNICAL SPECIFICATIONS

	GIMETTE 21	GIMETTE 28
Inside dimensions mm:	300x250x260	380x250x260
Approx. volume	20 dm ³	28 dm ³
Outside dimensions mm:	440x340x400	550x340x400
Weight:	13 kg	17 kg
Power consumption:	600 Watt	500 Watt
Operating Voltage:	230 V - 50/60 Hz	
Fuses:	2 x 4 A	
Timer:	manual from 0 to 120 minutes	
Thermostat:	automatic	
Thermometer:	high precision with capillary bulb	
2° Thermostat:	bi-metallic certified	
Green LED:	timer, orange LED: heating, red LED: alarm	

SPARE ACCESSORIES

STERILIZER	TRAY	GASKET	RESISTOR 600 W
Gimette 21	35610	35614	35639
Gimette 28	35611	35615	35639
Gimette 50*	-	35617	35639

*Needs 2 resistors