

SOAPED SPONGE, SOAPED AND NON SOAPED WASH GLOVE



36660

• **36660 SOAPED SPONGE** - box of 500 Disposable sponge in foam soaked with anallergenic, neutral PH soap. Ideal for patient cleaning and hygiene. Individually packed. Size: 7.5x12xh 2.5 cm



36663

• **36663 SOAPED WASH GLOVE** - box of 50 Disposable mitten to wash patient. Made in non woven and sponge, soaked with liquid soap with neutral PH and anallergenic. Ideal for patient cleaning and hygiene. Size: 15x25 cm



36661

• **36661 WASH GLOVE** - box of 1,000 Disposable wash glove to clean patients. Made in special airlaid, which allows to soap and rinse patient, and waterproof PE. Embossed surface allows a deep cleaning. Size: 16.5x22 cm. Made in Italy.

SKIN SOAP

• **36671 FRAGRANCE SKIN SOAP** - 1,000 ml

Delicate skin protective. Fragrance soothing bath and shower with moisturizer active agents, respects the lipid mantle and it is delicate on skin leaving it fresh and pleasantly soft. It is used for daily hygiene as a normal liquid soap, suitable for personal hygiene, bath, shower and shampoo. Suitable for frequent washing and very sensitive skins too, it protects and soothes sensitive skin of children and elderly people. Fast to be rinsed leaves a fresh cleanness sensation. Sold in box of 12 pcs, only 2 pieces come with dispenser pump.



• **36690 DISPENSER PUMP**

GERMELLA

• **36672 GERMELLA** - 1,000 ml

Germella is a skin-protecting liquid detergent combining the properties of three extraordinary medical plants: **Camomile, Sage and Aloe**. Camomile and Sage and have moisturizing and soothing qualities. Aloe Vera's main property is soothing. Moreover Aloe makes your skin soft and velvety. Thanks to the Camomile, Sage and Aloe extracts, Germella is a delicate detergent, which respects the natural defence of skin and mucous membranes. Thanks to its natural ingredients Germella is also particularly suitable for daily intimate care. Sold in box of 12 pcs, only 2 pieces come with dispenser pump.



BIBS WITH POCKET

Size: 38x28 cm + 10 cm (pocket)

36666

Easy fastening on patient by laces

Extremely soft and comfortable

Large pocket



• **36666 BIBS WITH POCKET** - box of 600 Disposable bibs made of 1 ply of highly absorbent wadding of cellulose and 1 ply of waterproof blue polyethylene. Box of 600 (6 bags of 100). Made in Italy.

HOT-COLD GEL

Individually packed in single box

34114



• **34114 HOT-COLD THERMO GEL** Size: 14x28 cm - box of 30 pcs. Special gel keeps temperature (hot or cold) for long time.

HOT WATER BOTTLES

Capacity 2,000 cc



28600

28601

• **28600 HOT WATER BOTTLE** - British standard - blue 100% made of rubber, latex free.
• **28601 HOT WATER BOTTLE 2 l** - red 25% made of rubber, latex free. Both are two side ribbed. Packed in plastic bag with multilingual user manual (GB, FR, IT, DE, ES, PT, GR).

ICE BAGS

Neck of bag: PP

Neck of bag: aluminium



28605

28606

• **28605 ICE BAG** - Ø 28 cm - large top cap O.D. 8.78 cm, I.D. 6 cm - light blue Made of polyester fabric + PVC, latex free. Single packed in carton box, with multilingual manual (GB, FR, IT, DE, ES, PT, GR).
• **28606 ICE BAG** - Ø 28 cm - small top cap O.D. 7.15 cm, I.D. 4.7 cm - blue Made of polyester fabric + PVC. Single packed in transparent bag with multilingual manual (GB, FR, IT, DE, ES, PT, GR)

ELECTRIC HEATING PAD

• **28655 ELECTRIC HEATING PAD** Suitable for many applications.
• **28656 ELECTRIC "SAND" HEATING PAD**

The use of purified marine sands make it possible to naturally transmit their beneficial warmth. Heating: position 1 up to 42°C in 30 min, position 2 up to 66°C (max 56°C first hour, up to 66°C following hours) External cover in cotton (94%) and nylon, internal bag in PVC Operating power: 230 V - 50/60 Hz - 75 W (28656) - 55 W (28655)



Sold in 4 languages box (GB, FR, IT, ES)

Size: 38x30 cm