

Sonotrax Vascular



Sonotrax Vascular offers a traditional pencil-style 4 or 8MHZ probe. It can locate both superficial and deep vessels faster and provides stronger signals.

The hand-held Doppler has been widely adopted as standard equipment for the assessment of peripheral vascular disease. The versatile, inexpensive and portable equipment allows the fast and simple but effective assessment of both arterial and venous disease. Besides, it helps physicians to diagnose some diseases, such as lower limb and foot ulceration in both diabetics and non-diabetics alike.

Sonotrax Vascular

Ultrasonic Pocket Doppler

Features

- Easy-to-use
- Ergonomic design
- Interchangeable and waterproof probe
- Rapid and effective

Sonotrax Vasular

Ultrasonic Pocket Doppler



specifications

Audio output power: 0.5W

Auto shut-off: 1 minute

Size: 32(Depth) x 85(Width) x 138(Height)mm

Weight: 290 ± 5g(including one battery)



Edan Instruments, Inc.

3F-B, Nanshan Medical Equipments Park, Nanhai Rd 1019#, Shekou, Nanshan, Shenzhen, 518067 P.R. China
Tel: 86-755-2689 8326 Fax: 86-755-2689 8330 www.edan.com.cn Email: info@edan.com.cn



All rights reserved.

Features and specifications are subject to change without notice

ENG-VA-vascular-V1.1-20090429