

## GYNAECOLOGICAL AND DELIVERY BEDS



Size: 180x55xh 80 cm  
Seat size: 55x55 cm  
Delivered in kit form.

27501

Multilingual instructions:  
GB, FR, IT, ES, PT, DE, GR.

Standard leg holders



360° rotation  
27503

Quality leg holders



360° rotation  
27504

MADE IN ITALY



Size: 180x55xh 80 cm.  
Delivered in kit form.

27510

- 27500 DELUXE GYN BED - quality leg holders
  - 27501 GYNAECOLOGY BED - standard leg holders
- Strong chrome-plated structure. Rectangular steel tube. Adjustable head and leg sections. Simulated leather, black colour, washable upholstery, self extinguishing. Load: 110 kg. Complete with s/s bowl Ø 31 cm and leg-holders.
- 26586 S/S BOWL Ø 31 cm - spare

- 27510 DELUXE DELIVERY BED - quality leg holders
- Same as 27500, but with three sections with independent controls and Trendelenburg position. Load: 110 kg.
- 27503 STANDARD LEG HOLDERS - couple - spare
  - 27504 QUALITY LEG HOLDERS - couple - spare
- Kit leg holders includes 2 leg-holders and arms.
- 27418 COUCH ROLL HOLDER

## 3-SECTION MULTIFUNCTIONAL AND GYNAECOLOGICAL BEDS



LOAD 150 KG

27505

Adjustable standard leg holder, turnable 360° (27508)



27509

Seat: 50x60 cm

Size: 181x60xh 83 cm  
Delivered in kit form.

MADE IN ITALY

- 27502 3-SECTION MULTIFUNCTIONAL BED - blue
  - 27505 3-SECTION GYNAECOLOGICAL BED - blue
- Strong and stylish structure in oval tube 50x25 mm, painted with epoxy powder, metallic grey colour. The three sections top has high density padding, covered with washable, fire proof, leatherette material Class 1 IM. Back-rest and leg-rest section are fully adjustable by self mechanical system. Gynaecological bed is supplied with 2 foam padded leg holders and stainless steel bowl Ø 32 cm. Maximum load: 150 kg - Weight: 36 kg.



Delivered in kit form

Seat: 62x60 cm

27502

Size: 185x60xh 85 cm

- 27508 STANDARD LEG HOLDER for 27505 - couple - spare
  - 27509 QUALITY LEG HOLDER for 27505 - couple - optional
- Kit leg holder, includes 2 leg-holders and arms.
- 27511 COUCH ROLL HOLDER for 27502-5
  - 27919 DESK with 2 drawers - 130x71xh 80 cm
- Strong structure in oval tube 40x20 painted with epoxy powder, metal grey colour. White work top with rubber edge and round corner.
- 27979 FOOT STOOL - 2 steps



27919



27979



27538

Height adjustable leg holder, turnable 360°

## 3-SECTION DELIVERY BED WITH TRENDELENBURG

- 27538 3-SECTION DELIVERY BED - with Trendelenburg
- structure made of epoxy painted oval steel tube
  - three sections with independent control and Trendelenburg position
  - padding made of high density polyurethane foam, fireproof Class 1 IM
  - covering made of washable, plastic material (leatherette) fireproof Class 1 IM
  - adjustable rubber feet

Size: 181x60xh 80 cm  
Carton size: 75x20xh 196 cm.  
Maximum load: 150 kg.  
Net weight: 43 kg.

Supplied with a pair of height-adjustable, padded, anatomic leg-holders, a couple of shoulder rests, handgrips, a paper roll holder and a stainless steel basin Ø 32 cm. Delivered in kit form. In conformity with MDR.

- 27509 QUALITY LEG HOLDER - couple - spare

MADE IN ITALY