

# 3M™ Littmann™ Lightweight II S.E. Stethoscope



Now with patented  
3M™ Littmann™  
Snap Tight Soft-  
Sealing Eartips for  
improved acoustics

**New features!**

Three new tubing  
colours give you a total  
of six colour choices

Ceil Blue	New
Seafoam Green	New
Lilac	New
Caribbean Blue	
Burgundy	
Black	

New chestpiece  
shape designed for  
easier positioning  
and assessment

Now with tunable  
diaphragm for easy  
low- and high-  
frequency monitoring

**A general-purpose  
stethoscope with  
an improved design**

# 3M™ Littmann™ Lightweight II S.E. Stethoscope




## Welcome to Littmann value.

This general purpose, entry-level stethoscope is an excellent, affordable choice for basic blood-pressure and limited physical assessment. Littmann quality shines through from the soft-sealing eartips down to the combination bell/diaphragm chestpiece that provides clear, reliable acoustic performance.

- New chestpiece shape designed for better orientation around blood-pressure cuffs and body contours
- Patented 3M™ Littmann™ Snap Tight Soft-Sealing Eartips provide a comfortable fit and a better seal to reduce ambient noise and improve acoustics
- Tunable diaphragm conveniently alternates between low- and high-frequency sounds without the need to turn over the chestpiece
- Combination bell/diaphragm chestpiece provides clear, reliable acoustic performance
- Lightweight design for easy portability
- Anatomically angled headset ensures proper fit for better acoustics
- Double-leaf binaural spring allows easy tension adjustment for comfortable fit
- Patient-friendly, nonchill rim and diaphragm
- Single-lumen tubing available in a variety of colours (including three new colours)
- Two-year warranty

### Tubing Colours



### Ordering Information

Catalog #	Length	Tube Colour
<b>Littmann Lightweight II S.E. Stethoscope</b>		
2450	28 in. (71.1 cm)	Black
2451	28 in. (71.1 cm)	Burgundy
2452	28 in. (71.1 cm)	Caribbean Blue
2453	28 in. (71.1 cm)	Lilac
2454	28 in. (71.1 cm)	Ceil Blue
2455	28 in. (71.1 cm)	Seafoam Green



Health Care  
3M Canada

Post Office Box 5757  
London, Ontario N6A 4T1  
1 800 364-3577 (8:00 a.m. - 6:00 p.m)  
www.3M.com/Littmann

For healthcare professional use only.  
Registered under Act 737.  
MDA Reg. No. : GMD74701541318A

3M and Littmann are trademarks of 3M. Used under license in Canada.  
© 3M 2004 All Rights Reserved  
0312-CP-18807-E