Tecnosix PR-TH

INTERMITTENT PNEUMATIC COMPRESSION

TECNOSIX PR-TH is an air compression system suitable to temporarily relieve muscle tension and increase blood circulation so as to promote the rapid recovery of the treated areas.



SUPPLIED WITH

- 1 control unit
- 2 legs
- 1 set of connecting pipes
- 1 adapter
- 1 carrying bag



TECHNICAL FEATURES:

- Chambers: 4
- Compression mode: 6
- Time range: 1 min - 99 min
- Adjustable pressure: 30 - 150 mmHg
- Rechargeable battery
- Adapter: yes (15V 6A)
- Input voltage: AC 100 -240V 50 - 60 Hz
- Dimension: cm 27 x 15 x 13
- Weight: 1,65 Kg (without battery)
- Sleeve garment: 420 denier

AVAILABLE ON REQUEST

- Sleeve Small: M165 cm
- Sleeve Medium: da 170 a 185 cm
- Sleeve Large: U185 cm
- Abdomen Sleeve

CE 1639

TECNOSIX PR-TH - medium size	**PRS0001**
TECNOSIX PR-TH - small size	**PRS0002**
TECNOSIX PR-TH - large size	**PRS0003**

