Continuous passive motion *for elbow*

XRI006 FISIOTEK LT-G

The Fisiotek LT-G is a mobilizer used for the passive rehabilitation of the elbow, by means of flexion-extension and pronosupination movements. The structure of the device provides great stability and precision of movement to ensure effective recovery. The programming of the range of movement is done electro-mechanically. The graphic scale on the central disc makes programming easy and intuitive. It is possible to adjust the height of the unit; the movement can be supervised by remote control with Start & Stop. Quality mechanical parts ensure the correctness of the movement over time.





Accessories set included

CE 0123



Elbow prono-supination exercise



Elbow flexion-extension exercise

Therapeutic indications

- Arthrotomy, arthroscopy
- Treatment following mobilization
- Surgical treatment of fractures
- Reconstructive surgery
- Endoprosthetic implants
- Operations on soft tissue

TECHNICAL DATA		
Elbow	Flexion-extension	0° ÷ 150°
	Pronation-supination	90° ÷ 0° ÷ 90°
Speed		min. 2°/sec. ÷ max. 4,5°/sec.
Range of movement		Adjustment by electromechanical limit switch
Power supply		100÷240 V AC - 50÷60Hz 400mA
Electric safety		Class II B Standard EN 60601-1
Electromagnetic compatibility		Group 1 Class B, Standard EN 60601-1-2
Classification according to EEC Directive 93/42		Class IIa
Weight		47 kg
Dimensions		75 x 75 x 103 h cm

Optional Accessories



02115 BATTERY CHARGER