Wide loop

MD CE

We offer both regular and wide loops. Wide loops are ideal for spreading forces more evenly on your body parts. They are high quality, cost effective and easy to use. The Moves Wide Loops are multi-layered and available in different resistances, thus satisfying all users.





SPECIFICATIONS

- Lay-flat-length: 30cm
- Width: 7,5cm

Colour	Resistance	Reference code 10-pack
	Tan | Extra light	01-100071
	Yellow | Light	01-100072
	Red | Medium	01-100073
	Green | Heavy	01-100074
	Blue | Extra heavy	01-100075
	Black | Special heavy	01-100076

PACKAGING

- 1 unit: 10 individually packed pcs of the same colour/resistance

---

Loop

MD CE

Moves Loops are extremely portable and unbelievably versatile products. They are very effective instruments for stretching and strengthening all muscle groups. The Moves Loops are multi-layered, which makes them thicker and more durable. The benefit of these thicker loops is that they won't roll up during exercising, which avoids the uncomfortable feeling from thinner alternatives.



SPECIFICATIONS

- Lay-flat-length: 30cm
- Width: 2,5cm
- Contains natural latex rubber

PACKAGING

- Polybag: 10 individual pcs of the same colour/resistance.
- Mesh bag: Set of 4 pcs of the most popular resistances

Colour	Resistance	Reference code 10-pack
	Tan | Extra light	01-100001
	Yellow | Light	01-100002
	Red | Medium	01-100003
	Green | Heavy	01-100004
	Blue | Extra heavy	01-100005
	Black | Special heavy	01-100006
	Set of 4 resistances	01-100011

	Total Kg / Elongation %					
	50%	**100%**	**150%**	**200%**	**250%**	**300%**
	1,9	2,8	3,4	4,0	4,7	5,3
	2,4	3,5	4,3	5,0	5,8	6,7
	3,1	4,5	5,6	6,6	7,6	8,7
	3,8	5,6	6,9,	8,1	9,3	10,7
	4,8	7,0	8,6	10,1	11,7	13,3
	6,0	8,7	10,7	12,6	14,6	16,7