

VT2

Professional floor scale

- Platform designed for weighing small and medium sized animals
- Large, easy to read LCD screen
- Handle and wheels built-in to enable transport
- Display can be wall mounted

Functions
Semi-automatic tare, net/gross weight, hold weight block, reset to zero, zero tracking, programmable auto-off, backlight activation, keyboard external calibration, low battery alert, overload signal

Structure
Painted metal tube structure with washable painted weighing plate (IRON) or with washable stainless steel weighing plate (INOX). 2m long power cable

Dimensions
W900 x D550 x H70 mm

Power supply
Rechargeable battery with 45h run time or power cable 15V-300mA

Packaging
W1000 x D600 x H100 mm
Weight: net 17,8 kg - gross 21 kg



VT2 is equipped with integrated wheels for an easier movement



Non-slip mat for extra stability and safety while weighing

Model	Code	Capacity	Division	Surface
VT2 IRON	00090V	150/300 kg	50/100 g	W900 x D550 mm
VT2 INOX	00090VA			

Included accessories: **power supply, rechargeable battery, RS232, adjustable feet, adjustable display support, 2 years warranty, non-slip mat**

Optional Services	Code
Calibration certificate w/ Accredia weights	ACOSIT