
## Ever-Trac Specifications:

| | |
|---|---|
| Traction force | The maximum traction force is 200 lb. (approximately 91kg)±4.5lb (2kg) |
| Force display | Actual force /preset force - automatic switching display |
| Treatment time | 1-99 minutes (±2%) |
| Safety system | Multiple protection alarms |
| Power supply(Volt) | 100-120V (60Hz), 220-240V (50Hz) |
| Electric safety degree abd complied standard | Class I/Type BF (EN60601-1,IEC 60601-1) |
| EMC complied standard | EN60601-1-2,IEC 60601-1-2 |
| CE Marking | CE 0434 |
| Weight | 30.8lb (14kg) |
| Dimensions | 31(L) x 28(W) x 26(H) cm |

## Variable Height Traction Couch Specifications:

| | |
|---|---|
| Minimum height | 470mm / 18.1" |
| Maximum height | 910mm / 35.8" |
| Cushion width | 700mm / 27.6" |
| Head section length | 410mm / 16.1" |
| Mid section length | 600mm / 26.8" |
| Leg section length | 740mm / 29.14" |
| Total length | 2240mm / 88.3" |
| Head tilt angles | +80° to-20° |
| Leg tilt angles | +87° to-0° |
| Lift mechanism | Electric / Hydraulic |
| Net weight | 113.9kg |


(1)


(2)


(3)


(4)


(5)

## Optional Extras

**Model#**

(1) **CA265** Head halter
(2) **CA260** Thoracic traction belt
(3) **CA255** Lumbar traction belt
(4) **CA270** Spreader bar
(5) **CA275** Counter Traction Straps

EVERYWAY

ISO 10993 SGS Approved


EVERYWAY MEDICAL INSTRUMENTS CO., LTD.
3F., No.5, Ln . 155, Sec. 3, Beishen Rd.,Shenkeng Dist., New Taipei City 22203,Taiwan
Telephone:+886 2 26620038 | Fax: + 8662 2664 5566
Web:www.everyway-medical.com | E-Mail:everyway@ms2.hinet.net

EVERYWAY

# Ever-Trac
# Lumbar & Cervical Traction

## Model:ET-800

**Continuous and intermittent traction, the ET-800 uses the latest technology combined with extremely simple operation and functional design. For both lumbar and cervial traction.**




CE 0434

# Everyway Variable Height Traction Couch

## Ever-Trac Features

**Model Number: ET-800**

- Crystal clear displays
- Static,Intermittent,Harmonic Traction
- Variable speed control
- Minimum and maximum force settings 2kgs /4.5lb - 91kgs/200lbs
- Patient safety switch with audible alarm
- Traction force compensation
- Audible alarm protection
- Continuous electronic monitoring of traction forces
- CE marked to the medical device directive 93/42/EEC

**Upholstery Colours**



Black Blue Beige Grey



Flexion Stool

## Contents of Standard Kit

(1) Traction unit x 1
(2) Power Cord x 1
(3) Patient Safety Switch x 1
(4) Fuses x 2
(5) Cord Extension x 1
(6) Quick Release Screws x 4
(7) Instruction Maunal x 1
(8) Service Manual x 1




(7)




(8)



(1)



(6)



(5)



(4)



(3)



(2)

## Standard Accessories

(1) **CA 251** Flexion Stool X 1
(2) **CA 254** Cervical Master Guide X 1
(3) **CA 285** Cervical In-sitting Attachment W/Cord Extension X 1

## Ordering Information

| Series | Model# | |
|---|---|---|
| S4 | EU200 | Electrical Variable height Traction Couch w/flexion stool(Angled or straight pull Option) |
| F4 | CA250 | Electrical Variable height traction couch w/flexion stool (Straight pull only) |
| F4 | CA255 | Hydraulic traction couch |

## Traction Couch Features

The Everyway Traction Couch is designed to offer the therapist multiple choices of traction therapy. The four piece top offers friction free movement for precise lumbar traction therapy. Cervical traction can be achieved in both lying and sitting positions. The cervical in sitting attachments, flexion stool are all supplied as standards.

Tel:+886 2 2662 0038 | Fax: +886 2 2664 5566 Email:everyway@ms2.hinet.net | Webmail.http://www.everyway-medical.com/