
Ozone-therapy





Vetozon compact: for the ozone therapy in the veterinary sector

Excellent operational comfort and high safety – these are the characteristics of systems by HERRMANN Apparatebau. Vetozon compact has been developed according to our strict quality guidelines and can be used for almost the entire range of modern ozone treatment. The integrated colour touch screen allows for an uncomplicated handling.

Highlights of this system:

- Can be used for all types of normobaric ozone therapy
- High flexibility thanks to the possiblility to draw up ozone into a syringe
- High compatibility thanks to Luer Lock connections
- Multilingual operational guidance (German, English, Spanish, Italian, Turkish)
- Colour touch screen for a more comfortable handling

Technical data:

Devicename:	Vetozon compact
Power supply system:	230V~/50-60 Hz or 115V~/50-60 Hz
Input:	240 VA
Concentration ozone:	0 and 5-80 µg/ml
Flow through:	0,8 l/min
Dimensions (BxHxT):	424 x 185 x 388 mm
Weight:	13,1 kg

☐ 230V Ord. no. 25002
☐ 115V Ord. no. 25003

Vetozon compact

HERRMANN Apparatebau GmbH • Dieselstraße 8 • 63839 Kleinwallstadt • Germany
+49 (0) 6022-6581-3 • +49 (0) 6022-6581-59
www.h-a-b.de • info@h-a-b.de