iHealth Bluetooth glucose monitor - see page 3

# GLUCOSE MONITORS FOR DAILY DIABETES TESTING

Advanced technology glucose monitors. Available as unit only or complete kit.
Kit includes 10 strips, code chip, control solution, 10 sterile lancets, lancet device, clear cap (for testing on forearm or palm), carrying case, user manual.
Main features:
- accurate result in just 5 seconds
- only 0.5 µl sample size
- automatic calibration with code chip
- plasma calibration results
- advanced data management with 7, 14, 30 days average
- affordable for frequent monitoring
- flexibility to test on the forearm, palm or fingertip
- certified ISO 15197:2013

CE 0123

- BIOSENSOR TECHNOLOGY
- 5 SECOND TEST TIME
- ONLY 0.5 µL BLOOD SAMPLE

*Easy sample application with capillary action*




23891

*Easy calibration with chip PLASMA calibration results*

| TECHNICAL SPECIFICATIONS | |
|---|---|
| Measurement Range | 20 to 600 mg/dL or 1.1 to 33.3 mmol/L |
| Sample Sites - Size | Finger tips, Forearm or palm - 0.5 µL |
| Test Time | Approximately 5 seconds |
| Power Source | One 2032 or equivalent 3.0 V coin cell battery |
| Battery Life | 12 months or approximately 1,000 tests |
| Memory | Up to 300 records with time and date |

| mg/dL | | | mmol/L | ON CALL® PLUS II GLUCOSE MONITOR |
|---|---|---|---|---|
| Italian Greek | English French | Spanish Portuguese | English French | |
| 23890 | 23891 | 23892 | 23881 | Glucose monitor Plus-meter only |
| 23895 | 23896 | 23897 | 23886 | Glucose monitor Plus-kit |

| GIMA code | Accessories and consumables for GLUCOSE MONITOR ON CALL® PLUS AND ON CALL® PLUS II |
|---|---|
| 23910 | Glucose Strips - box of 25 strips |
| 23912 | Glucose Strips - box of 50 strips |
| 23913 | Glucose Strips - box of 100 strips |
| 23914 | Control Solution |
| 23919 | Auto lancet device |
| 23917 | Lancets 28G for 23919 - sterile - box of 25 pcs |
| 23916 | Lancets 28G for 23919 - sterile - box of 100 pcs |




23896

# ABBOTT GLUCOSE STRIPS




Abbott

| ABBOTT FREESTYLE OPTIUM BLOOD GLUCOSE TEST STRIPS | | |
|---|---|---|
| GIMA code | Compatibility | Minimum order |
| 23942 | FreeStyle Optium Neo, FreeStyle Libre, FreeStyle Optium, Optium Xceed, Optium, Optium Easy, Boots | box of 25 |
| 23943 | | box of 50 |
| **ABBOTT FREESTYLE LITE BLOOD GLUCOSE TEST STRIPS** | | |
| 23944 | FreeStyle Lite, FreeStyle Freedom Lite, FreeStyle InsuLinx, FreeStyle Navigator II | box of 25 |
| 23945 | | box of 50 |

# LANCETS AND LANCET DEVICE

**• 23919 AUTO LANCET DEVICE**
Reusable auto lancet plastic device, adjustable tip, allows precise configuration of lancet penetration.



23916




23917

**• 23917 GIMA LANCETS - 28G - sterile - box of 25 pcs**
**• 23916 GIMA LANCETS - 28G - sterile - box of 100 pcs**
Lancets are fine needles that can be used with code 23919 and most of lancet devices in the market to draw a blood sample for glucose testing. Multilanguage box: GB, FR, IT, DE, ES, PT, GR, PL, Arabic.

23916-17



# SURGILANCE SAFETY AUTOMATIC LANCETS

**SURGILANCE SAFETY AUTOMATIC LANCETS - NEEDLES**
Multi purpose capillary blood sampling devices for all testing requirements.
The one-step safety lancet is easier to use, no arming is required, and safer, once lancet is used it is rendered inoperable.
High-speed delivery and penetration method minimizes patient pain and operator error.
Latex-free, hypo-allergenic, disposable, sterile.

| SURGILANCE™ SAFETY LANCETS - NEEDLES | | | | | | | | | |
|---|---|---|---|---|---|---|---|---|---|
| GIMA code | Colour | Depth of penetration | Entire needle Gauge | Needle tip Gauge | Blood Flow | Finger Stick | Glucose | Hematocrit | Minimum order |
| 24513 | Yellow | 1.0 mm | G21 | G26 | 5-10 µl | • | • | | box of 100 |
| 24514 | Grey | 1.8 mm | G21 | G23 | 10-20 µl | • | • | | box of 100 |
| 24515 | Orange | 2.2 mm | G21 | G22 | 20-40 µl | • | • | • | box of 100 |
| 24517 | Pink | 2.8 mm | G21 | G21 | 40-60 µl | • | • | • | box of 100 |


*Multilanguage manual on www.medipurpose.com*




CE 1008

*Multilanguage box: GB, FR, IT, ES, DE, PT, NL, SE, DK, RU, GR.*

24515



INSULINE PEN NEEDLES SEE PAGE 52