WUNDER PROFESSIONAL LABORATORY SCALE WJ600 M



Product Code	31030
Unit of sale	1 pc
Minimum order	1
Type	
Class	

UK-REP	No
CH-REP	No

EAN/UPC	8023279310306

Description

WUNDER PROFESSIONAL LABORATORY SCALE WJ600 M

Compact, reliable and portable laboratory scale, designed for rapid and accurate precision weighing. ABS structure with a waterproof keyboard and washable stainless-steel weighing plate. Large 20mm backlit LCD display with 5 digits.

Functions: net/gross weight, 10x high-resolution mode, subtractive tare, piece counting, weight %, 7 units (g, ct, oz, ozt, dwt, t, gn), external digital calibration, overload alarm, auto-zero and programmable auto-off.

Supplied with power supply, rechargeable battery, RS232, level, adjustable feet, windshield.

Made in Italy.

Technical Specifications

Capacity - sensitivity: 600 g - d=0,01 g e=0,1 g
Weighing surface: Ø 115 mm
Power supply: rechargeable battery with 40h run time or 12V power cable
Size - weight: 200x225xh 80 mm - 1.3 kg
Packing size: 320x320xh 120 mm