FERROX®
Physiotherapy Quality Solutions

CEZANNE JR 5

Physiotherapy table with five sections.



Medical
Device

Dimensions
197 x 64 x 50 h92 cm

Adjustable height
from cm 50 min to 92 cm max

Safe working load
160 Kg

Maximum capacity
200 Kg

Unimpeded passage
for patient lifting

⚠
The photo depicts the product complete with optional accessories

FEATURES

- Electric motor with a thrust force of 8000 N
- Safety control key
- Push-button control panel
- Class 1IM fire-retardant padding and upholstery