
ENDOSCOPY PRODUCTS

Our endoscopy portfolio combines patented rod lens technology and premium German optics to deliver crystal-clear images and an exceptionally bright field of view. Crafted with stainless-steel bodies, scratch-resistant sapphire lenses, and intuitive directional markings, these instruments offer durability. All endoscopes include STORZ, WOLF, OLYMPUS, and ACM connectors. Additionally, the entire range is fully compatible with both low-temperature plasma and EO sterilization.


LAPAROSCOPES



- **71995 LAPAROSCOPE Ø 10x325 mm - 0° - autoclavable**
- **71996 LAPAROSCOPE Ø 10x325 mm - 30° - autoclavable**

- Micro-nano ultra-smooth surface processing technology, scratch and corrosion resistance
- 70° field angle to widen the visual scope


SINUSCOPES



- **72014 SINUSCOPE Ø 4x175 mm - 0° - non autoclavable**
- **72016 SINUSCOPE Ø 4x175 mm - 70° - non autoclavable**

Made of stainless steel with good corrosion resistance.
60° field of view.

OTOSCOPES



- **72001 OTOSCOPE Ø 2.7x110 mm - 0° - non autoclavable**
- **72002 OTOSCOPE Ø 3x110 mm - 30° - non autoclavable**

Made of stainless steel with good corrosion resistance.
50° field of view.

4K LAPAROSCOPE



- **71991 4K FLUORESCENCE LAPAROSCOPE Ø10x340 mm - 30° - autoclavable**

- to transmit white light and near infrared exciting light spectrum
- the aspheric optical technology greatly reduces the image edge distortion, improving image resolution
- 4K definition, providing the richest source optical information
- high temp. and pressure sterilization shorten process time
- 75° large viewing angle, large depth of field

LARYNGOSCOPE



- **72022 LARYNGOSCOPE Ø 8x184 mm - 70° - non autoclavable**

Made of stainless steel with good corrosion resistance.
45° field of view.

RESECTOSCOPES



- **72031 RESECTOSCOPE Ø 2.9x302 mm - 30° - non autoclavable**
- **72032 RESECTOSCOPE Ø 4x302 mm - 30° - autoclavable**

- the stainless-steel water valve radically eliminates the water valve maintenance procedure and weakness of easy damage
- when the inner and outer sheaths remain still, the manipulator, endoscope and surgical electrode can rotate 360° for continuous perfusion
- 60° field of view

NEUROLOGICAL REFLEX HAMMERS

GIMA code	HAMMERS - BLACK EDITION	Colour	Handle	Pin	Brush	Weight	Length	Multilingual box (all) and instructions
31266	Buck - black edition	black ■	stainless steel	yes	yes	100 g	18 cm	GB, FR, IT, ES, PT, DE, Arabic
31267	Dejerine - black edition	black ■	stainless steel	yes	yes	150 g	22 cm	Instructions available via QR code: NL, DK, SE, NO,
31268	Babinsky - black edition	black ■	stainless steel	yes	no	160 g	24 cm	FI, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR
GIMA code	**LUXAMED HAMMERS**	**Colour**	**Handle**	**Pin**	**Brush**	**Weight**	**Length**	**Multilingual label and instructions**
31248	Buck	white □	aluminium/brass	yes	no	90 g	18 cm	GB, FR, IT, ES, PT, DE, PL
31250	Buck	magenta ■	aluminium/brass	yes	no	90 g	18 cm	GB, FR, IT, ES, PT, DE, PL



