e-xam

Diagnostic penlight



Photobiologically validated, high-performance penlights for medical examination of the throat and pupils. Available with LED or Xenon illumination.

- Photobiologically validated for pupil illumination in accordance with EN 62471
- Pre-focused, uniform light field
- Bright, natural white light from high-performance LED
- LED lifespan of at least 20,000 hours
- Low maintenance costs thanks to long LED and battery life
- ABS casing with metal clip; disinfectable; easy to use with ON/OFF switch
- Supplied with blade holder, batteries, and nylon case
- Designed for durability; supplied with a 2-year warranty, including LED



Technical data

e-xam **LED 2.5 V**	5000 - 5500 K
	180,000 Lux max.
	field of light: approx. 11 cm diameter (distance 15 cm)
e-xam **Xenon 2.5 V**	3500 - 4000 K
	60,000 Lux max.
	field of light: approx. 14 cm diameter (distance 15 cm)
Batteries	2 x Alkaline AA, 1.5 V

ORDERING INFORMATION

5130-01	e-xam, XL 2.5 V, black, with tongue blade holder
5131-01	e-xam, LED 2.5 V, black, with tongue blade holder

906140 Rev. A 2026-01-26

Please visit **riester.de** for more information

Rudolf Riester GmbH P.O. Box 35 | Bruckstraße 31 | 72417 Jungingen | Germany
Tel.: (+49) 7477-9270-0 | Fax.: (+49) 7477-9270-70 | info@riester.de | www.riester.de