ri-pen®

Diagnostic pupil penlight



Sustainable penlight for medical diagnosis with LED illumination.

- Photobiologicaly validated for lighting of pupils according to EN 62471
- New design with pupil gauge from 2 mm - 9 mm and a ruler with 8 cm [3 in]
- UV free light through LED illumination in an aluminium casing with practical metal clip
- Supplied as single units or in packs of 6 (assorted colors or 6 per color)
- Simple ON/OFF function with push button on the end cap of the penlight
- Comes with 2 AAA batteries at time of purchase

Technical data

ri-pen® LED	3 V
	5300 K
	40,000 Lux max.
Field of view	Approximately 40 mm diameter at a working distance of 15 cm
Batteries	2 x Alkaline AAA, 1.5 V

ORDERING INFORMATION

5075-526	**ri-pen®** black, pack of 1 pc.	**5075**	**ri-pen®** black, pack of 6 pcs.
5074-526	**ri-pen®** silver, pack of 1 pc.	**5074**	**ri-pen®** silver, pack of 6 pcs.
5077-526	**ri-pen®** red, pack of 1 pc.	**5077**	**ri-pen®** red, pack of 6 pcs.
5078-526	**ri-pen®** purple, pack of 1 pc.	**5078**	**ri-pen®** purple, pack of 6 pcs.
5071-526	**ri-pen®** blue, pack of 1 pc.	**5071**	**ri-pen®** blue, pack of 6 pcs.
5072-526	**ri-pen®** green, pack of 1 pc.	**5072**	**ri-pen®** green, pack of 6 pcs.
		5070	**ri-pen®** assorted colors, pack of 6 pcs.

906130 Rev. A 2026-01-26

Please visit **riester.de** for more information

Rudolf Riester GmbH P.O. Box 35 | Bruckstraße 31 | 72417 Jungingen | Germany
Tel.: (+49) 7477-9270-0 | Fax.: (+49) 7477-9270-70 | info@riester.de | www.riester.de