SC Latodis com SRL
Poiana Campina, jud. Prahova
Str. Centru nr.636,bl.21,parter
Tel. (+40) 744 511 725
E-mail: office@latodis-med.com

SuperField®






VSFNC

PRIMARY APPLICATION

Wide Field Small Pupil Pan Retinal Examination

- The 'Super 90D' – same magnification as the 90D with a wider field of view enabling both posterior pole and pan retinal examinations
- Provides dynamic, high resolution viewing to the periphery
- Combines a wide field of view with a comfortable working distance and magnification
- Can be used on small pupils and patients who do not accommodate dilation

95°/116°	0.76x	1.32x
FIELD OF VIEW	IMAGE MAG	LASER SPOT MAG



Lentile de diagnostic SuperField® Volk, tehnologie avansata, recomandate în diagnosticul oftalmologic modern. Lentila diagnostic Volk® optica de calitate superioara, ofera imagini detaliate in cadrul examinarii campului periferic si vizualizare dinamica si rapida a retinei.

Lentila SuperField®, camp vizual mare, inel lentila de 30 mm, prindere confortabila si manipulare usoara, combinatie perfecta a campului vizual larg si a distantei de lucru confortabila asigura utilizarea perfecta la pacientii (pupile mici) care nu accepta dilatarea. Lentilele Volk® reprezinta alegerea ideala pentru profesionistii din domeniul oftalmologiei în evaluarea si diagnosticul afectiunilor oculare.

Lentile diagnostic:SuperField® Volk

Camp vizual: 95° / 116°

Marire imagini: 0.76x

Marire spot laser: 1,32x

Distanta de lucru: 7 mm

Culori disponibile: negru, la cerere - albastru, auriu, verde, violet, rosu, argintiu







