SC Latodis com SRL
Poiana Campina, jud. Prahova
Str. Centru nr.636,bl.21,parter
Tel. (+40) 744 511 725
E-mail: office@latodis-med.com

G-3 Gonio



Flange: **VG3**
No Flange: **VG3NF** (shown)

Gonio Mini
No flange: **VG3MININF** (shown)


Available in mini version for pediatric and patients with small orbits

PRIMARY APPLICATION

3-Mirror, All-Glass Design for Anterior, Peripheral, and Equatorial Viewing (Goldmann-style Lens)

- 60° mirror provides a view of the iridocorneal angle to visualize the trabecular meshwork
- 66° mirror provides a retinal image from the equator to the ora serrata
- 76° mirror provides a view of the mid-peripheral/far-peripheral retina
- Central lens enables clear viewing of the posterior pole
- Available in two formats: flanged (can be used for laser, however, not compatible with SLT) and no flanged (recommended for routine gonioscopy without indentation)

60°/66°/76°	1.06x	0.94x
MIRROR ANGLES	IMAGE MAG	LASER SPOT MAG



Lentila diagnostic Volk® Gonio G-3, alegerea ideala pentru profesionistii din domeniul oftalmologiei în evaluarea si diagnosticul oftalmologic al afectiunilor oculare. Lentile de diagnostic Volk® Gonio design cu trei oglinzi din sticla asigura o vizualizarea clara si detaliata a unghiului iridocornean, a retinei central periferice si a retinei medii periferice/periferice. Cu o calitate imbatabila a opticii Volk lentila în stil Goldman este instrumentul perfect pentru specialistii în glaucom.

Lentila diagnostic Volk® Gonio G-3 cu flansa pentru vizualizare anterioara, periferica si ecuatoriala. Oglinda de 60˚ ofera o vedere a unghiului iridocornean, oglinda de 66˚ ofera o imagine retiniana central-periferica, oglinda de 76˚ ofera o vedere a retinei medii periferice/periferice. Lentila centrala ofera o vedere clara a polului posterior

Lentile diagnostic: Volk® Gonio G-3

Camp vizual: 60°/66°/76°

Marire imagini: 1,06x

Marire spot laser: 0,94x

Distanta de lucru: 15 mm

Culori disponibile: albastru







