
SC Latodis com SRL
Poiana Campina, jud. Prahova
Str. Centru nr.636,bl.21,parter
Tel. (+40) 744 511 725
E-mail: office@latodis-med.com



Lentila de diagnostic Super Pupil® XL, tehnologie avansata, recomandata în diagnosticul oftalmologic pentru examinarea retiniana pupilei mici de 2-3 mm. Lentila diagnostic Volk® optica de calitate superioara, ofera imagini detaliate, camp vizual mare, combinatia perfecta a campului vizual larg si a distantei de lucru confortabila asigura utilizarea perfecta la pacienti. Lentilele Volk® reprezinta alegerea ideala pentru profesionistii din domeniul oftalmologiei în evaluarea si diagnosticul afectiunilor oculare.

Lentile diagnostic: Volk Super Pupil® XL

Camp vizual: 103° / 124°

Marire imagini: 0.45x

Marire spot laser: 2.20x

Distanta de lucru: 4 mm

Culori disponibile: negru, la cerere - albastru, auriu, verde, violet, rosu, argintiu







