Wenus | Gynaecology examination chair

Mars | Urological examination chair





Wenus
A chair designed for gynaecology examination and treatment.

Mars
A chair designed for urological examination and urodynamic testing.

One chair, so many advantages!

Electric motors for easy adjustment of seat and backrest, including Trendelenburg position. Motors controlled with foot switch or hand-operated remote controller. Braked castors. Paper roll holder behind the backrest. Special bowl holder that folds under the seat. Optional panel turns the chair into a versatile operating table.



Main features:

- versatile
- ergonomic
- robust
- mobile
- smooth, precise adjustment
- high-quality finish

Technical data:

Length (excl. operating table panel) [cm]:	180 (135)
Width (excl. clamps) [cm]:	90 (78)
Seat width [cm]:	60
Backrest width [cm]:	60
Seat adjustment (chair mode) [cm]:	47-100
Trendelenburg position adjustment [°]:	20 (+/-2)
Motor voltage [VDC]:	24
Work cycle of motors:	2 min. work / 18 min. pause
Allowable load [kg]:	220
Weight (exc. accessories) [kg]:	~100
Power [VAC]:	230
Upholstery colour:	sea green (default), other colours available

Upholstery colours:

00 ecru	01 beige	02 blue	03 navy blue
04 terracotta	05 brown	06 grey	07 graphite
08 black	09 red	10 light green	11 honey
12 sea green	13 orange	14 lime	15 light blue

Accessories for Wenus chair:





Leg rest with a GIN clamp mount



Bowl without a drain

Shared frame of Mars and Wenus chairs



Operating table panel




Arm rest with an URO clamp mount for drip infusions



Common accessories

Accessories for Mars chair:





Leg rests with an URO clamp mount



Hand rests



Swivel leg rests



Bowl
(with and without a drain)



Flowmeter funnel holder



Wenus/Mars: options and accessories:

View	Catalogue number	Description
	01-1505-0300-0000 01-1520-2000-0004	Gynaecology examination chair WENUS + leg rests with GIN clamp mount + paper roll holder. No controllers included. Purchased separately.
	01-1505-0300-0000 01-1520-1900-0003 01-1520-2000-0005	Urological examination chair MARS + leg rests with URO clamp mount + hand rests + paper roll holder. No controllers included. Purchased separately.
	01-1520-0900-0002	Arm rest (left or right) with URO clamp mount for drip infusions (1 pc).
	01-1520-0600-0005	Operating table panel. Turns both chairs into typical operating tables.
	01-1520-0600-0008	Swivel leg rests. Swivel leg rests for MARS for urodynamic testing. Set (2 pcs).
	01-1520-0300-0003	Urological (MARS) bowl with drain.
	01-1520-0300-0004	Urological or gynaecology (MARS/WENUS) bowl without drain.
	01-1520-0300-0005	Flowmeter funnel holder.
	01-1520-2600-001	Infusion hanger.
	01-1520-1500-0004	Foot-operated controller.
	01-1520-1800-0002	Hand-operated remote controller.

MEDEN INMED
PROFESSIONAL MEDICAL ENGINEERING

MEDEN-INMED Sp. z o.o., ul. Wenedów 2, 75-847 Koszalin
Phone: +48 94 344 90 78, fax: +48 94 347 10 41
km@meden.com.pl, tgorecki@meden.com.pl
www.meden.com.pl