AR-EL TIBBI ÜRÜNLER

2087 GYNECOLOGICAL OBSTETRIC DELIVERY TABLE

with 4 Motors Electrically adjustable height – backrest – trendelenburg – Knee Crutches
with 3 Motors Electrically adjustable height – backrest – trendelenburg manual Knee Crutches





Techinical Details:
Seat Height **:** 640mm - 940 mm
Table Lenght **:** 2010mm
Table Width **:** 640mm
Seat Section **:** 640mm x 360mm
Back Section **:** 640mm x 1000mm
Footrest Section : 640mm x 600mm
Curved Leg Section **:** Min space:850mm
Lifting Capacity **:** 200 Kg
Motor **:** 4 or 3 pieces
Command Module **:** Foot switch





Functional Details:
Back Section **:** 0° ~ +50°
Trendelenburg **:** -5° ~ +25°
Foot Extension **:** Detachable
Curved Legrest **:** Adjustable upwards-downwards & left- right and Detachable

Standart Accessories :
Foot Extension (1 pc), Curved Legrest (2 pcs), Knee crutches (2 pcs), Handle (2 pcs), Waste tray (1 pc), Paper roll holder (1 pc), Accessory hanger (1 pc).

Optional Accessories :
Foot Step (2 pcs)

Packing Details:
Net weight**:** 139 Kg
Gross weight **:**169 Kg
Packing Dimensions **:** 152 x 77 x 99 cm
Packing Type **:** Full protect wooden case.

Warranty : 2 year

