FERROX®
Physiotherapy Quality Solutions

CHAGALL NEO 6

Six-section electric osteopathy table Chagall Neo 6 is suitable for global postural re-education as well as active and passive treatments based on myofascial release techniques and postural gymnastics. It is designed with variable height and backrest tilt assisted by electric motor.



Medical
Device

Dimensions
196 x 65 x 49-94h cm

Adjustable height
da cm 49 cm min a 94 cm max

Safe working load
250

Maximum capacity
320

Unimpeded passage
for patient lifting

⚠
Version presented with perimeter command and wheel options

FEATURES

"Camel hump" and Trendelenburg positions

Electric motor with a thrust force of 8000 N for height

Electric motor with a thrust force of 60000 N for the backrest

Push-button control panel

Safety control key

Footrest tilt assisted by two gas pistons

Class IM1 fire-retardant padding and upholstery, washable