FERROX®
Physiotherapy Quality Solutions

CEZANNE PRO 3

Osteopathy table with three segments and very compact foam



Medical
Device

Dimensions
190 x 54 x 50 h92 cm

Adjustable height
da cm 52 cm min a 96 cm max

Safe working load
150 Kg

Maximum capacity
190 Kg

Unimpeded passage
for patient lifting

⚠
The photo depicts the product complete with optional accessories

FEATURES

- Three-section table with gas piston-assisted forward and backward tilting of the backrest
- "Camel hump" and Trendelenburg positions
- Electric motor with a thrust force of 4000N for the "camel hump" function
- Electric motor with a thrust force of 8000 N for elevation
- Safety control key
- Class IM1 fire-retardant padding and upholstery