FERROX®
Physiotherapy Quality Solutions

CEZANNE PRO 5

Four-segment osteopathy table with electric height and kyphosis adjustment. Rapid lifting time of 43 seconds.



Medical
Device

Dimensions
193 x 64 x 52-96h cm

Adjustable height
da cm 52 cm min a 96 cm max

Safe working load
150

Maximum capacity
190

Unimpeded passage
for patient lifting

⚠
Version presented with perimeter command and wheel options

FEATURES

- 100° pivoting armrests on the headrest
- "Camel hump" and Trendelenburg position assisted by electric motor
- Electric motor with a thrust force of 4000N for the "camel hump" function
- Electric motor with a thrust force of 8000 N for elevation
- Push-button control panel
- Safety control key
- Class IM1 fire-retardant padding and upholstery