
FERROX®
Physiotherapy Quality Solutions

CEZANNE PRO 7

Seven-segment electric osteopathy table with adjustable height and kyphosis-assisted function. Rapid lifting time of 43 seconds.



Medical
Device

Dimensions
194 x 64 x 52-96h cm

Adjustable height
da cm 52 cm min a 96 cm max

Safe working load
150

Maximum capacity
190

Unimpeded passage
for patient lifting

⚠
Version presented with perimeter command and wheel options

FEATURES

180° pivoting armrests on the headrest

"Camel hump" and Trendelenburg position assisted by electric motor

Electric motor with a thrust force of 4000N for the "camel hump" function

Electric motor with a thrust force of 8000 N for elevation

Foldable pectoral armrests on the central section

Push-button control panel

Safety control key

Class IM1 fire-retardant padding and upholstery