# EXAMINATION AND TREATMENT TABLE - HIGH LOAD


- **27622 STANDARD TREATMENT TABLE - any colour****
- **27623 STANDARD TREATMENT TABLE - blue***
- **27624 STANDARD TREATMENT TABLE - cream***
- **27626 LARGE TREATMENT TABLE - any colour****
- **27627 LARGE TREATMENT TABLE - blue***

Professional examination and treatment table with one movable section and manual adjustment.
Sturdy detachable metal structure painted with scratch-proof, non corrosive epoxy powder.
The section trim facilitates reaching the patient.
The head section has an extended mouth-nose hole.
Load 180 kg (standard), 200 kg (large).
Made in Italy

- **27629 COUCH ROLL HOLDER**
- **27630 MOUTH-NOSE PLUG - any colour****
- **27631 MOUTH-NOSE PLUG - blue**
- **27620 MOUTH-NOSE PLUG - cream**

** usually in stock  ** available in 30 days*

AVAILABLE IN 6 COLORS




*1 Black  4 Green*
*2 Blue  5 Cream*
*3 Grey  6 Apricot*




**27622**

*High-density foam padding. Superior cover finish made of high-quality fireproof simulated leather. Vinyl is resistant to disinfectants, urine, and blood.*





**27627**

WIDTH 80 cm

MAX LOAD - 200 kg

# WOODEN COUCHES

**WOODEN COUCHES**
Comfortable couches, ideal for physical therapy, cardiology and wellness centres. Superior quality wood. Back rest adjustable in 13 positions. Resilient foam padding covered with washable weave. Padding in beige colour.
Size: 190x70xh 70 cm.



**27416**




**27417+27409**
*Mouth nose hole plug included*

*Maximum load: 150 kg*

- **27416 WOODEN COUCH - beech**
- **27417 WOODEN COUCH - with face hole - beech**
- **27443 DANTE WOODEN TROLLEY - beech**
- **27409 COUCH ROLL HOLDER for 27416 and 27417**

Couch roll holder included with 27450 and 27451. Delivered in kit form. Easy assembly, tool provided. Made in Italy.

*Backrest adjustable in 18 positions*

*Mouth nose hole plug included*

**27451**



**27452**
*Size 50x40xh 77 cm*



**27450**

*Maximum load: 200 kg*

- **27450 FIRENZE WOODEN COUCH - wenge**
- **27451 FIRENZE WOODEN COUCH - with face hole - cherry wood**
- **27452 DANTE WOODEN TROLLEY - cherry wood**

# SUITCASE COUCH

- **27407 SUITCASE COUCH**

Chromed steel and folding legs.
Upholstered covering in washable material.
Resin padded. Available in black colour.
Weight: 19 kg - Max load: 95 kg
Supplied with a velvet bag for easy transport.
Made in Italy.



**27407**

*Couch size:*
*185x64xh 70 cm (open)*
*95x65x16 cm (folded)*

# ALUMINIUM SUITCASE COUCH

- **27406 ALUMINIUM SUITCASE COUCH with hole (22x12)**

Light, foldable, two sections, suitcase couch with aluminium frame. Central double legs for improved safety and block system on the hinges. Head rest with mouth/nose hole.
Upholstered and covered with simulated leather.
Maximum load: 110 kg.
Weight: 13 kg. Made in Italy.



*Size: 62x182xh 75 cm*

**27406**