FERROX®
Physiotherapy Quality Solutions

PICASSO XL

Two-section physiotherapy electric table with forward tilting of the backrest only, and electrically assisted height adjustment. Recommended for therapy sessions, massages and medical exams.



Medical
Device

Dimensions
194 x 64 x 50 h92 cm

Adjustable height
da cm 50 cm min a 92 cm max

Safe working load
160 Kg

Maximum capacity
200 Kg

Unimpeded passage
for patient lifting

⚠
Version presented with perimeter command and wheel options

FEATURES

- Two-section table with with 84 cm forward tilting of the backrest only, assisted by gas piston.
- Electric motor with a thrust force of 8000 N
- Push-button control panel
- Safety control key
- Class 1IM fire-retardant padding and upholstery