BOBATH JR 200x120


Treatment table with electrically adjustable height control for Bobath and Vojta therapy
Versions shown with optional castors



Detail: reinforced, independent backrest frame

1






(1) **BOBATH Jr - Art. 240006 TE -**
As described for Art. 240006 but with electrically assisted backrest

- Hand switch control
- Medical device

See detailed features page 57




OPTIONS



240006 TER


2





(2) **BOBATH Jr - Art. 240006**
- Electric motor with a thrust force of **8000 N**
- Foot pedal control
- Safety control key
- Fire-retardant padding and covering class 1IM
- Electric height adjustment range 48 Δ 90 cm
- Safe working load 160 kg
- Dimensions: 200 x 120 x 48 Δ 90 cm
- Medical device

See detailed features page 56



OPTIONS



240006 R

BOBATH JR 200x100


Treatment table with electrically adjustable height control for Bobath and Vojta therapy
Versions shown with optional castors




1

(1) **BOBATH Jr - Art. 240007**
- Electric motor with a thrust force of **8000 N**
- **Manual** control switch
- Safety control key
- Fire-retardant padding and covering class 1IM
- Electric height adjustment range 48 Δ 90 cm
- Safe working load 160 kg
- Dimensions: 200 x 100 x 48 Δ 90 cm
- Medical device

See detailed features page 57





OPTIONS



240007 R





2

(2) **BOBATH Jr - Art. 240007 TE -**
As described for Art. 240007 but with electrically assisted backrest

- Hand switch control
- Medical device

See detailed features page 57





OPTIONS



240007 TER