# CLINEA NEBULIZER - piston




| | |
|---|---|
| Product Code | 28058 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | II A |

| | |
|---|---|
| UK-REP | No |
| CH-REP | Yes |

| | |
|---|---|
| EAN/UPC | 6970384250077 |
| GTIN | 06970384250077 |

## Description

CLINEA NEBULIZER - piston

High-end, low-noise compressor nebulizer designed to provide efficient and quiet aerosol therapy.
Ideal for professional healthcare settings.
Featuring innovative brushless DC motor technology, this device combines power, reliability and long-lasting performance to deliver a comfortable and effective treatment experience.
Equipped with an adjustable pressure function, the nebulizer allows users to adjust the airflow intensity according to their specific needs or medical recommendations.

Supplied with adult mask, child mask, bulb, inhalation mouthpiece, inhalation tube, air filter, adapter and multilingual manual GB, IT

## Technical Specifications

Operating mode: continuous use
Nebulizing flow: ≥4 L/min
MMAD: <5 µm
Noise level: ≤55 dB
Expected pump life: 1000 hours
FPF: ≥60%
Nebulizing Rate: ≥0.2 ml/min
Max capacity: 10 ml
Operating time: > 4 hours
Size: 168.5x116x99 mm
Voltage: 100~240V 50~60Hz, 700mA Max
Input: DC 12V 1.5 A