

# GIMA ENT HEAD MIRRORS FOR BATTERY OR AC/DC TRANSFORMER USE

HEAD MIRRORS WITH TRANSFORMER



Head mirror Lux-100 with plug in transformer




HEAD MIRRORS WITH BATTERY







LUX-100 AVAILABLE WITH LED LIGHT

Quality head mirrors Ø 55 or 100 mm to perfectly illuminate the examination area in every situation.
Long life 6 V lamp provides bright and natural light, quality silver coated mirror allows optimum light reflection.
Mirror and light beam can be easily and precisely adjusted.

| STANDARD ACCESSORIES | LUX AC/DC TRANSFORMER | LUX BATTERY* | |
|---|---|---|---|
| | 31110-31111 | 31107-31109 | 31112 |
| Mirror with head band | • | • | • |
| Long life bulb | • | • | • |
| Head band mounted battery case** for 2 x 3V CR-123A | | • | |
| Plug in transformer* 100/240 V - 50/60 Hz with cable | • | | • |
| Pocket battery case | | | • |
| Sturdy plastic case | | | • |
| Cordura bag | • | • | |
| Extension cable 80 cm and on/off switch 30 cm cable | | | • |

*Supplied with interchangeable UK, US, EU plugs **Battery not included

| GIMA code | LUX ENT MIRRORS |
|---|---|
| 31108 | LUX-55 battery** head mirror - Ø 55 mm - 6 V |
| 31112 | LUX-55 battery*** head mirror - Ø 55 mm - 6 V |
| 31109 | LUX-100 battery** head mirror - Ø 100 mm - 6 V |
| 31107 | LUX-100 battery** LED head mirror - Ø 100 mm - 6 V |
| 31110 | LUX-55 AC/DC transformer head mirror - Ø 55 mm - 6 V |
| 31111 | LUX-100 AC/DC transformer head mirror - Ø 100 mm - 6 V |
| 31114 | Mirror for LUX-55 - spare |
| 31115 | Mirror for LUX-100 - spare |
| 31116 | Bulb for LUX mirrors - spare |
| 31119 | Carrying bag for LUX mirrors - black |
| 31117 | Battery case for LUX mirror 31112/3 - spare |
| 31124 | Connecting cable 30 cm for 31107-9 - spare |
| 31118 | Plug in transformer* 100-220 V - 50/60 Hz with cable - spare for 31100-1 - optional for 31107-9 |


*Supplied with interchangeable UK, US, EU plugs **Battery mounted on head mirror band
***Pocket battery case connected to head mirror

# HEAD LOUPE AND MAGNIFYING GLASSES




**• 31170 HEAD LOUPE two magnifying systems (2.2X and 3.3X)**
The system is mounted on an adjustable webbing head band and provided with two separate lenses. The first lens provides a 2.2X magnification, moving down the second lens, we reach a total 3.3X magnification. The superior quality of this lens permits high precision examinations. It can be used in combination with any kind of glasses. Working distance: 17 cm with 2.2X lens, 9 cm with 3.3 lens. Made in Japan.

**• 30861 DANAE MAGNIFYING GLASSES**
*Light, only 83 g*




*Compatible with TETI protective shield (see page 16, 198)*

- optical system with 3 magnifying glasses (2X, 2.5X, 3X)
- can fit over your glasses
- adjustable to any face shape

**• 30848 SMART MAGNIFYING GLASSES**
Spectacle frame with magnifying glasses useful to physician whenever it is necessary to have a magnified view of the visual field and be able to use both hands. Lightweight, it folds up, and can be worn over spectacles. 2.5X zoom lens which can be tilted to achieve optimum focus.




# ZIEGLER MIRROR




**• 31142 ZIEGLER MIRROR - Ø 90 mm**
Concave back-silvered mirror, 90 mm diameter with 13 mm diameter hole, on adjustable band. Special mount permits replacement of glass mirror in event of breakage.
Made in Germany.

# 1 W LED HEADLIGHTS

**• 30858 LED HEADLIGHT**
Run time: High mode: 3 hours - Low mode: 22 hours
Luminous intensity: 120 lumens (high) or 30 lumens (low)
30 feet drop test, super bright CREE cool white LED
MCU control push switch (High/Low/Off) 60° pivoting head.
Passed extreme temperature shock test, IPX4 water resistant
10,000 switch cycle life
Lighting leakage design when forget to turn off switch. 3AAA batteries included.

**• 32200 "AAA" ALKALINE BATTERIES - blister of 4 pcs - spare**




*Compact and lightweight headlight only 65 g*