DIAGNOSTIC INSTRUMENTS
Ophthalmology, ENT & Dermatology



ri-scope® L

High performance for life

The ri-scope® L - Premium handheld
LED Diagnostic instruments



Riester

ri-scope® L Otoscope

Premium otoscopes with glass optics and LED / Li-ion technology



ri-scope® L2

The economic model with 2.5 V / 3.5 V LED fiber optic illumination, lasting up to 100,000 hours (without SpecEjec specula ejection device)

- Bi-directional swivelling optical glass with 3x magnification, optional operation lens with 4x magnification available for purchase
- Metal fitting for **e**-scope® / **ri**-mini® specula, compatible with specula by other well-known manufacturers
- With integrated fiber optic throat illuminator
- Suitable for pneumatic otoscopy (supplied without pneumatic bulb)
- Large selection of power sources: rechargeable Li-ion batteries, disposable alkaline batteries or lithium batteries; desktop charging base or for use with ri-former diagnostic stations.


ri-scope® L3

The high-end version with 2.5 V / 3.5 V LED fiber optic illumination, lasting up to 100,000 hours with SpecEjec specula ejection device and operation lens.

- Fiber optics for optimal transmission of the light
- Bi-directional swivelling optical glass with 3x magnification
- Includes operation lens with 4x magnification and glass mini optic ease the insertion of instruments
- With integrated fiber optic throat illuminator
- Easy-to-operate SpecEjec specula ejection device at the rear side of the otoscope to reduce cross contamination
- Fitting for **ri**-scope® specula and enlarged conduit for instruments
- Suitable for pneumatic otoscopy (supplied without pneumatic bulb)
- Large selection of power sources: rechargeable Li-ion batteries, disposable alkaline batteries or lithium batteries; desktop charging base or for use with **ri**-former diagnostic stations.



Otoscope
Head with operation lens



The SpecEjec device enables touchless ejection of used ear specula for hygienic and simplified handling.

ri-scope® L Otoscope Head

10582*	**ri**-scope® L2 otoscope head, LED 2.5 V	**10583***	**ri**-scope® L3 otoscope head, LED 2.5 V
10565	**ri**-scope® L2 otoscope head, LED 3.5 V	**10567**	**ri**-scope® L3 otoscope head, LED 3.5 V
10565-301	**ri**-scope® L2 otoscope head, LED 3.5 V with anti-theft security	**10567-301**	**ri**-scope® L3 otoscope head, LED 3.5 V with anti-theft security

*Compatible with alkaline batteries only and the **ri**-scope® L handle (**10670** – sold separately). Not suitable for rechargeable batteries or use with **ri**-former®.

ri-scope® L Ophthalmoscope

Premium-Class coaxial ophthalmoscope with renowned high-performance optics and LED / Li-ion technology




ri-scope® L2


All the important basic functions with separately engageable red-free, blue and polarization filters for each aperture.

- Dioptre wheel with 29 corrective lenses **Plus** 1 - 10, 12, 15, 20, 40 **Minus** 1 - 10, 15, 20, 25, 30, 35
- Easy-to-operate aperture hand-wheel with slit, fixation star, semi-circle and small, medium and large circles
- Includes filter wheel engageable for all apertures with symbol display, red-free filter, blue filter and polarization filter
- Optimized high-performance optics with aspherical condenser lens
- Parallel beam path
- Dust-proof
- Large selection of power sources: rechargeable Li-ion batteries, disposable alkaline batteries or lithium batteries; desktop charging base or for use with ri-former diagnostic stations.

Apertures L2



Filters



ri-scope® L3


The maximized model with additional focussing device and 89 single-step dioptre values.

- Corrective lenses for 89 dioptre values **Plus** 1 - 45 in single steps **Minus** 1 - 44 in single steps
- Easy-to-operate aperture hand-wheel with grid, slit, fixation star, semi-circle and small, medium and large circles
- Includes filter wheel engageable for all apertures with symbol display, red-free filter, blue filter and polarization filter
- Focusing device with zero retention force for the quick correction of ametropias
- Dust-proof
- Large selection of power sources: rechargeable Li-ion batteries, disposable alkaline batteries or lithium batteries; desktop charging base or for use with ri-former diagnostic stations.




Ophtalmoscope Head

Apertures L3



Filters



ri-scope® L Otoscope Head

10597*	**ri-scope®** L2 ophthalmoscope head, LED 2.5 V	**10598***	**ri-scope®** L3 ophthalmoscope head, LED 2.5 V
10571-203	**ri-scope®** L2 ophthalmoscope head, LED 3.5 V	**10573-203**	**ri-scope®** L3 ophthalmoscope head, LED 3.5 V
10595-301	**ri-scope®** L2 ophthalmoscope head, LED 3.5 V with anti-theft security	**10596-301**	**ri-scope®** L3 ophthalmoscope head, LED 3.5 V with anti-theft security

*Compatible with alkaline batteries only and the **ri-scope®** L handle (**10670** – sold separately). Not suitable for rechargeable batteries or use with **ri-former®**.

DIAGNOSTIC INSTRUMENT HEAD

Riester diagnostic instruments are engineered for precision, reliability and comfort. Featuring advanced 3.5 V LED technology with a lifespan of up to 100,000 hours, they deliver cool, daylight-like illumination through high-quality fiberoptic (F.O.) lightwave conductors. Each instrument is designed with ergonomic and durable casings — ranging from impact-resistant plastics to chrome-plated metal — and offers quick, secure bayonet fittings and easy lamp replacement. Whether examining skin, throat, or nasal passages, these tools ensure optimal visibility and a seamless clinical experience.



EliteVue Macro Otoscope

- More than 5x magnification
- Complete view of the tympanic membrane
- Reflection-optimised, high-performance glass optics
- Wheel to adjust focal depth (+5 to -5)
- Easy-to-operate SpecEjec specula ejection device at the rear side of the otoscope to reduce cross contamination

10510*	**Elite**Vue head, LED 2.5 V
10512	**Elite**Vue head, LED 3.5 V
10512-301	**Elite**Vue head, LED 3.5 V with anti-theft security



ri-derma dermatoscope

- 10x magnification with focusing achromatic lens
- Stable black chrome-plated metal casing, dust-proof
- Practical eyeglass protection

10577	**ri-**scope® **ri-**derma dermatoscope head, LED 3.5 V
10577-301	**ri-**scope® **ri-**derma dermatoscope head, LED 3.5 V with anti-theft security



F.O. tongue blade holder

- Slide attachment for secure fixing of wooden tongue depressors
- Sturdy, lightweight ABS plastic casing

10574	**ri-**scope® F.O. tongue blade holder, LED 3.5 V
10574-301	**ri-**scope® F.O. tongue blade holder, LED 3.5 V with anti-theft security



F.O. nasal speculum

- Swiveling and removable lens with 2.5x magnification
- Controlled expansion via screw adjuster; fast expansion by pressing down
- Impact-resistant glass-fiber reinforced plastic casing

10575	**ri-**scope® F.O. nasal speculum, LED 3.5 V
10575-301	**ri-**scope® F.O. nasal speculum, LED 3.5 V with anti-theft security



F.O. light conductor

- Optimal illumination of throat and laryngeal areas
- Chrome-plated metal construction for durability

10576	**ri-**scope® F.O. light conductor, LED 3.5 V
10576-301	**ri-**scope® F.O. light conductor, LED 3.5 V with anti-theft security

*Compatible with alkaline batteries only and the **ri-**scope® L handle (**10670** – sold separately). Not suitable for rechargeable batteries or use with **ri-**former®.

MODULAR DIAGNOSTIC SYSTEM

Diagnostic station with flexible, modular design

ri-former® – a modular diagnostic station offering flexibility for hospitals and medical practices. The diagnostic modules can be combined individually to suit the practitioner's specialty and location. The **ri-**former® diagnostic station exhibits notable advantages in terms of its versatile configurations.

With its centralised power supply, the **ri-**former® provides the necessary reliability to meet the demanding and fast-paced environment found in hospitals and medical practices.

Blood Pressure Modules

- Fully automated, clinical-grade blood pressure monitor [**RBP-**100]
- Manual, large-scale, clinical-grade aneroid sphygmomanometer [big ben®]

Thermometry Modules

- Infrared tympanic thermometer [**ri-**thermo® tym**PRO**]
- Non-contact, clinical-grade infrared thermometer [**ri-**thermo® sensio**PRO**]
- Predictive oral / axillary / rectal thermometer [**RPT-**100]

Instrument Heads

- **Riester** diagnostic **ri-**former® systems support up to three handles, each compatible with six interchangeable heads: otoscope, ophthalmoscope, dermatoscope, tongue blade holder, nasal speculum or light conductor.

Accessories

- **ri-**spec integrated specula dispenser – for ear and nasal specula, choose from five funnels for five different specula sizes. This allows the use of a new speculum for each examination and is easy to refill and disinfect.
- Extension module – for one additional handle



For **ri-**former® instrument heads with anti-theft device option, please see pages 2 – 4.
For all other options, please see separate **ri-**former® brochure.

DIAGNOSTIC INSTRUMENT SETS

ri-scope® L otoscope Set

- ri-scope® L otoscope head
- Type C battery handle (batteries not included)
- A tube containing 10 disposable ear specula in assorted sizes ranging from 2 mm to 5 mm
- Operating lens with 4x magnification (included only in the L3 set)
- Deluxe hard case with safety catch



3703	**ri-scope®** L2 otoscope, LED 2.5 V / Type C battery handle for 2 Type C alkaline batteries or 1 **ri-accu®**
3704	**ri-scope®** L3 otoscope, LED 2.5 V / Type C battery handle for 2 Type C alkaline batteries or 1 **ri-accu®**
3703-550	**ri-scope®** L2 otoscope, LED 3.5 V / Type C battery handle for 1 **ri-accu®** L
3704-550	**ri-scope®** L3 otoscope, LED 3.5 V / Type C battery handle for 1 **ri-accu®** L

ri-scope® L ophthalmoscope Set

- ri-scope® L ophthalmoscope head
- Type C battery handle (batteries not included)
- Separate magnifying glass with 5x magnification
- Deluxe hard case with safety catch



3725-203	**ri-scope®** L2 ophthalmoscope, LED 2.5 V / Type C battery handle for 2 Type C alkaline batteries or 1 **ri-accu®**
3726-203	**ri-scope®** L3 ophthalmoscope, LED 2.5 V / Type C battery handle for 2 Type C alkaline batteries or 1 **ri-accu®**
3811-203	**ri-scope®** L2 ophthalmoscope, LED 3.5 V / Type C battery handle for 1 **ri-accu®** L
3812-203	**ri-scope®** L3 ophthalmoscope, LED 3.5 V / Type C battery handle for 1 **ri-accu®** L

ri-scope® L otoscope / ophthalmoscope Set

- ri-scope® L otoscope and ophthalmoscope head
- Type C battery handle (batteries not included)
- A tube containing 10 disposable ear specula in assorted sizes ranging from 2 mm to 5 mm
- Operating lens with 4x magnification (included only in the L3 set)
- Deluxe hard case with safety catch



3746.013	**ri-scope®** L2 otoscope, LED 2.5 V / L2 ophthalmoscope, LED 2.5 V / Type C battery handle for 2 Type C alkaline batteries or 1 **ri-accu®**
3746-203	**ri-scope®** L3 otoscope, LED 2.5 V / L2 ophthalmoscope, LED 2.5 V / Type C battery handle for 2 Type C alkaline batteries or 1 **ri-accu®**
3817-204	**ri-scope®** L2 otoscope, LED 3.5 V / L3 ophthalmoscope, LED 3.5 V / Type C battery handle for 1 **ri-accu®** L
3746-204	**ri-scope®** L3 otoscope, LED 3.5 V / L2 ophthalmoscope, LED 3.5 V / 2 Type C battery handles with 2 **ri-accu®** L and **ri-charger®** L

ENT **Kit**

A premium diagnostic set combining the **ri-**scope® L2 otoscope and ophthalmoscope, both equipped with high-performance LED 3.5 V illumination for precise and reliable examinations.



- **ri-**scope® L2 otoscope and ophthalmoscope (LED 3.5 V)
- USB Type-C battery handle with rechargeable Li-ion battery
- USB-C to USB-A charging cable, charging adapter and an interchangeable plug adapter
- A tube containing 10 disposable ear specula in assorted sizes ranging from 2 mm to 5 mm
- Durable zip case for secure storage and portability

3832	**ri-**scope® L2 otoscope, LED 3.5 V / L2 ophthalmoscope, LED 3.5 V / USB Type C battery handle with Li-ion battery, US plug adapter
3833	**ri-**scope® L2 otoscope, LED 3.5 V / L2 ophthalmoscope, LED 3.5 V / USB Type C battery handle with Li-ion battery, EU plug adapter
3834	**ri-**scope® L2 otoscope, LED 3.5 V / L2 ophthalmoscope, LED 3.5 V / USB Type C battery handle with Li-ion battery, UK plug adapter
3835	**ri-**scope® L2 otoscope, LED 3.5 V / L2 ophthalmoscope, LED 3.5 V / USB Type C battery handle with Li-ion battery, AU plug adapter

EliteVue **otoscope / ri-**scope® L **ophthalmoscope Set**



- **Elite**Vue otoscope head and **ri-**scope® L2 ophthalmoscope head
- Type C battery handle (batteries included or excluded depending on set)
- A tube containing 10 disposable ear specula in assorted sizes ranging from 2 mm to 5 mm
- Deluxe hard case with safety catch

2210-204	**Elite**Vue otoscope / **ri-**scope® L2 ophthalmoscope, LED 2.5 V / Type C battery handle for 2 Type C alkaline batteries or 1 **ri-**accu®
2210-203	**Elite**Vue otoscope / **ri-**scope® L2 ophthalmoscope, LED 3.5 V / Type C battery handle for 1 **ri-**accu® L
2212-203	**Elite**Vue otoscope / **ri-**scope® L2 ophthalmoscope, LED 3.5 V / 2 Type C battery handles / 2 **ri-**accu® L / **ri-**charger® L

ri-scope® **E.N.T.**



- **ri-**scope® L otoscope head, **ri-**scope® F.O. fiber optic nasal speculum and **ri-**scope® F.O. fiber optic tongue blade holder
- Type C battery handle (batteries not included)
- A tube containing 10 disposable ear specula in assorted sizes ranging from 2 mm to 5 mm
- Operation lens with 4x magnification
- Deluxe hard case with safety catch

3760	**ri-**scope® L3 otoscope, LED 2.5 V / **ri-**scope® F.O. nasal speculum / **ri-**scope® F.O. tongue blade holder / Type C battery handle for 2 Type C alkaline batteries or 1 **ri-**accu®
3760-550	**ri-**scope® L3 otoscope, LED 3.5 V / **ri-**scope® F.O. nasal speculum / **ri-**scope® F.O. tongue blade holder / Type C battery handle for 1 **ri-**accu® L

ri-scope® L perfect



- ri-scope® L otoscope head, ophthalmoscope head and ri-scope® F.O. fiber optic tongue blade holder
- Type C battery handle (batteries not included)
- A tube containing 10 disposable ear specula in assorted sizes ranging from 2 mm to 5 mm
- Operation lens with 4x magnification
- Deluxe hard case with safety catch

3754	**ri-scope®** L3 otoscope, LED 2.5 V / L2 ophthalmoscope, LED 2.5 V / **ri-scope®** F.O. tongue blade holder / Type C battery handle for 2 Type C alkaline batteries or 1 **ri-accu®**
3819-203	**ri-scope®** L3 otoscope, LED 3.5 V / L2 ophthalmoscope, LED 3.5 V / **ri-scope®** F.O. tongue blade holder / Type C battery handle for 1 **ri-accu®** L

ri-scope® L praktikant



- ri-scope® L otoscope head, ri-scope® F.O. fiber optic nasal speculum, ri-scope® F.O. fiber optic tongue blade holder and ri-scope® F.O. fiber optic light conductor
- Type C battery handle (batteries not included)
- A tube containing 10 disposable ear specula in assorted sizes ranging from 2 mm to 5 mm
- Operation lens with 4x magnification
- Two laryngeal mirrors No. 3 and No. 4
- Deluxe hard case with safety catch

3766	**ri-scope®** L3 otoscope, LED 2.5 V / **ri-scope®** F.O. nasal speculum / **ri-scope®** F.O. tongue blade holder / **ri-scope®** F.O. light conductor / Type C battery handle for 2 Type C alkaline batteries or 1 **ri-accu®**
3766-550	**ri-scope®** L3 otoscope, LED 3.5 V / **ri-scope®** F.O. nasal speculum / **ri-scope®** F.O. tongue blade holder / **ri-scope®** F.O. light conductor / Type C battery handle for 1 **ri-accu®** L

ri-scope® L de luxe



- ri-scope® L otoscope head, ophthalmoscope head, ri-scope® F.O. fiber optic nasal speculum, ri-scope® F.O. fiber optic tongue blade holder and ri-scope® F.O. fiber optic light conductor
- Type C battery handle (batteries not included)
- A tube containing 10 disposable ear specula in assorted sizes ranging from 2 mm to 5 mm
- Two laryngeal mirrors No. 3 and No. 4
- Deluxe hard case with safety catch

3772	**ri-scope®** L3 otoscope, LED 2.5 V / L2 ophthalmoscope, LED 2.5 V / **ri-scope®** F.O. nasal speculum / **ri-scope®** F.O. tongue blade holder / **ri-scope®** F.O. light conductor / Type C battery handle for 2 Type C alkaline batteries or 1 **ri-accu®**
3820-203	**ri-scope®** L3 otoscope, LED 3.5 V / L2 ophthalmoscope, LED 3.5 V / **ri-scope®** F.O. nasal speculum / **ri-scope®** F.O. tongue blade holder / **ri-scope®** F.O. light conductor / Type C battery handle for 1 **ri-accu®** L

ri-scope® ri-derma dermatoscope



- ri-scope® L ri-derma dermatoscope head
- Type C battery handle (batteries not included)
- Deluxe hard case with safety catch

3777	**ri-scope® ri-derma** dermatoscope LED 2.5 V / Type C battery handle for 2 C alkaline batteries or 1 **ri-accu®**
3777-550	**ri-scope® ri-derma** dermatoscope LED 3.5 V / Type C battery handle for 1 **ri-accu®** L

HANDLE

Type C Battery Handle



- Equipped with rheotronic® electronic light adjustment technology – simply rotate the switch ring clockwise or counterclockwise for smooth, precise brightness control to suit clinician preference
- Solid metal construction with black chrome plating for durability and premium finish
- Practical bayonet fitting for quick and secure instrument head changes
- Compatible with 3.5 V and 2.5 V **ri-scope**® instrument heads
- Works with 2 C alkaline batteries or 1 **ri-accu**® / **ri-accu**® L rechargeable battery
- Simple charging via **plug-in charger** (with **ri-accu**® L **[10694]**) or **charging base** (**ri-charger**® L with **ri-accu**® L **[10691]** or **ri-charger**® with **ri-accu**® **[10681]**)
- Easy battery replacement at the base of the handle
- Handle diameter: 28 mm
- Supplied without batteries

10670	Type C battery handle with rheotronic®, chrome-plated, for 2 Type C alkaline batteries or 1 **ri-accu**® / **ri-accu**® L

CHARGING OPTIONS

Charging bases



- **ri-charger**® L (white model): For all rechargeable handles with Li-ion batteries, 3.5 V
- **ri-charger**® (black model): For all rechargeable handles with NiMH batteries, 2.5 V and 3.5 V
- Available as table-top or wall-mounted unit
- Includes a practical storage compartment for two replacement lamps and ear specula
- Charging status displayed on device
  - Yellow: charging in progress
  - Green: charging completed

10705	**ri-charger**® L 3.5 V, wall support with EU adapter	**10700**	**ri-charger**® 2.5 V / 3.5 V, wall support with EU adapter
10705.001	**ri-charger**® L 3.5 V, wall support with UK adapter		
10705.002	**ri-charger**® L 3.5 V, wall support with AU adapter	**10700.006**	**ri-charger**® 2.5 V / 3.5 V, wall support with AU adapter
10706	**ri-charger**® L 3.5 V, wall support with US adapter	**10701**	**ri-charger**® 2.5 V / 3.5 V, wall support with US adapter
2212-203	**Elite**Vue otoscope / **ri-scope**® L2 ophthalmoscope, LED 3.5 V / 2 Type C battery handles / 2 **ri-accu**® L / **ri-charger**® L		
3746-204	**ri-scope**® L set with LED L3 otoscope/LED L2 ophthalmoscope, 2 C-handles and ri-charger® L		

Accessories

10713	Ear specula compartment cover for **ri-charger**® L, blue
10710	Wall support **ri-charger**® / **ri-charger**® L, with screws and plugs

Plug-in charger



- Operate with 3.5 V **ri-accu**® L **[10694]** rechargeable Li-ion battery
- Up to 9 hours operating time (LED) when fully charged
- Latest charging technology – less than 3 hours charging time
- Wide-range plug-in charger (120 – 240 V) with regional plug options available on request

10707	Plug-in charger for rechargeable Li-ion battery 3.5 V **ri-accu**® L **[10694]** for Type C battery handle, w/o battery, EU plug
10708	Plug-in charger for rechargeable Li-ion battery 3.5 V **ri-accu**® L **[10694]** for Type C battery handle, with battery, EU plug
10708.002	Plug-in charger for rechargeable Li-ion battery 3.5 V **ri-accu**® L **[10694]** for Type C battery handle, with battery, US plug

USB-C Handle



- Higher capacity – up to 10 hours operating time
- Battery status indicator on handle (bi-colour LED shows: fully charged, charging, or charge required)
- Universal USB Type-C charging technology
- Wide-range plug-in charger (120 – 240 V) with regional plug options available on request
- Fully charged in approximately 3.5 hours
- Charging status
  - Green: Battery fully charged
  - Flashes green: Battery charging
  - Orange: Battery critically low – charge required

10704	USB Type C battery handle, 3.5 V with rechargeable Li-ion battery, charging cable, charging adapter and EU plug adapter
10704.001	USB Type C battery handle, 3.5 V with rechargeable Li-ion battery, charging cable, charging adapter and AU plug adapter
10704.002	USB Type C battery handle, 3.5 V with rechargeable Li-ion battery, charging cable, charging adapter and UK plug adapter
10704.003	USB Type C battery handle, 3.5 V with rechargeable Li-ion battery, charging cable, charging adapter and US plug adapter

BATTERIES

10691	**ri-accu**® L rechargeable Li-ion battery Type C, 3.5 V for Type C battery handle and **ri-charger**® L
10694	**ri-accu**® L rechargeable Li-ion battery Type C, 3.5 V for Type C battery handle and plug-in charger
10681	**ri-accu**® rechargeable NiMH battery Type C, 2.5 V for Type C battery handle
10699	**ri-accu**® USB rechargeable Li-ion battery Type C, 3.5 V for Type C battery handle
10686	Type C alkaline battery 1.5 V, pack of 2 pcs.

LAMPS



LED 2.5 V

- For **ri-scope**® otoscope L2 / L3
- Pack contains 1 unit

10626	LED 2.5 V lamp, 1 pc.



LED 2.5 V

- For **ri-scope**® ophthalmoscope L2 / L3
- Pack contains 1 unit

10624	LED 2.5 V lamp, 1 pc.



LED 3.5 V

- For **ri-scope**® otoscope L2 / L3, Light conductor, nasal speculum and tongue blade holder
- Pack contains 1 unit

10625	LED 3.5 V lamp, 1 pc.



LED 3.5 V

- For **ri-scope**® ophthalmoscope L2 / L3
- Pack contains 1 unit

10627	LED 3.5 V lamp, 1 pc.

SPECULA

Reusable specula for L2

- Available in sizes: 2 mm, 2.5 mm, 3 mm, 4 mm and 5 mm
- Made of black plastic, each speculum clearly labelled by size
- Supplied in a polybag



10775	Reusable ear specula Ø 2.0 mm, pack of 10 pcs.	**10779**	Reusable ear specula Ø 2.5 mm, pack of 10 pcs.
10783	Reusable ear specula Ø 3.0 mm, pack of 10 pcs.	**10789**	Reusable ear specula Ø 4.0 mm, pack of 10 pcs.
10795	Reusable ear specula Ø 5.0 mm, pack of 10 pcs.		

Disposable specula for L2



- For single use – a hygienic option to prevent cross-contamination
- Made of black plastic, each speculum clearly labelled by size
- Supplied in a polybag

10772-532	Disposable ear specula Ø 2.0 mm, pack of 100 pcs.	**10773-532**	Disposable ear specula Ø 2.0 mm, pack of 500 pcs.
10774-532	Disposable ear specula Ø 2.0 mm, pack of 1000 pcs.	**10772-531**	Disposable ear specula Ø 2.5 mm, pack of 100 pcs.
10773-531	Disposable ear specula Ø 2.5 mm, pack of 500 pcs.	**10774-531**	Disposable ear specula Ø 2.5 mm, pack of 1000 pcs.
10772-533	Disposable ear specula Ø 3.0 mm, pack of 100 pcs.	**10773-533**	Disposable ear specula Ø 3.0 mm, pack of 500 pcs.
10774-533	Disposable ear specula Ø 3.0 mm, pack of 1000 pcs.	**10772-534**	Disposable ear specula Ø 4.0 mm, pack of 100 pcs.
10773-534	Disposable ear specula Ø 4.0 mm, pack of 500 pcs.	**10774-534**	Disposable ear specula Ø 4.0 mm, pack of 1000 pcs.
10772-535	Disposable ear specula Ø 5.0 mm, pack of 100 pcs.	**10773-535**	Disposable ear specula Ø 5.0 mm, pack of 500 pcs.
10774-535	Disposable ear specula Ø 5.0 mm, pack of 1000 pcs.		

14061-532	Disposable specula Ø 2.0 mm, pack of 100 pcs.	**14062-532**	Disposable specula Ø 2.0 mm, pack of 500 pcs.
14063-532	Disposable specula Ø 2.0 mm, pack of 1000 pcs.	**14061-531**	Disposable specula Ø 2.5 mm, pack of 100 pcs.
14062-531	Disposable specula Ø 2.5 mm, pack of 500 pcs.	**14063-531**	Disposable specula Ø 2.5 mm, pack of 1000 pcs.
14061-533	Disposable specula Ø 3.0 mm, pack of 100 pcs.	**14062-533**	Disposable specula Ø 3.0 mm, pack of 500 pcs.
14063-533	Disposable specula Ø 3.0 mm, pack of 1000 pcs.	**14061-534**	Disposable specula Ø 4.0 mm, pack of 100 pcs.
14062-534	Disposable specula Ø 4.0 mm, pack of 500 pcs.	**14063-534**	Disposable specula Ø 4.0 mm, pack of 1000 pcs.
14061-535	Disposable specula Ø 5.0 mm, pack of 100 pcs.	**14062-535**	Disposable specula Ø 5.0 mm, pack of 500 pcs.
14063-535	Disposable specula Ø 5.0 mm, pack of 1000 pcs.		

Reusable specula for L3



- Specially developed for use with the SpecEjec option, absolutely reliable handling
- Available in sizes: 2 mm, 3 mm, 4 mm, 5 mm and 9 mm
- The 9 mm speculum is suitable for use as a nasal speculum
- Made of black plastic, each speculum clearly labelled by size
- Supplied in a polybag

10800-532	Reusable ear specula Ø 2.0 mm, pack of 10 pcs.	**10800-533**	Reusable ear specula Ø 3.0 mm, pack of 10 pcs.
10800-534	Reusable ear specula Ø 4.0 mm, pack of 10 pcs.	**10800-535**	Reusable ear specula Ø 5.0 mm, pack of 10 pcs.
10800-539	Reusable nasal specula Ø 9.0 mm, pack of 10 pcs.		

Disposable specula for L3



- For single use – a hygienic option to prevent cross-contamination
- Made of black plastic, each speculum clearly labelled by size
- Supplied in a polybag

10801-532	Disposable ear specula Ø 2.0 mm, pack of 100 pcs.	**10802-532**	Disposable ear specula Ø 2.0 mm, pack of 500 pcs.
10803-532	Disposable ear specula Ø 2.0 mm, pack of 1000 pcs.	**10801-533**	Disposable ear specula Ø 3.0 mm, pack of 100 pcs.
10802-533	Disposable ear specula Ø 3.0 mm, pack of 500 pcs.	**10803-533**	Disposable ear specula Ø 3.0 mm, pack of 1000 pcs.
10801-534	Disposable ear specula Ø 4,0 mm, pack of 100 pcs.	**10802-534**	Disposable ear specula Ø 4.0 mm, pack of 500 pcs.
10803-534	Disposable ear specula Ø 4.0 mm, pack of 1000 pcs.	**10801-535**	Disposable ear specula Ø 5.0 mm, pack of 100 pcs.
10802-535	Disposable ear specula Ø 5.0 mm, pack of 500 pcs.	**10803-535**	Disposable ear specula Ø 5.0 mm, pack of 1000 pcs.
10801-539	Disposable nasal specula Ø 9.0 mm, pack of 100 pcs.	**10802-539**	Disposable nasal specula Ø 9.0 mm, pack of 500 pcs.
10803-539	Disposable nasal specula Ø 9.0 mm, pack of 1000 pcs.		

ACCESSORIES



Operation lens

- Optical glass with a small inner diameter for instrument insertion
- 4x magnification
- With pin connector
- For use with L2 and L3 otoscopes

11449	Operation lens for **ri-scope®** L otoscope



Laryngeal mirrors

- For **ri-scope®** F.O. light conductor
- Two laryngeal mirrors No. 3 (Ø 20 mm / 0.79 in) and No. 4 (Ø 22 mm / 0.87 in) for illumination and mirroring of laryngeal area
- Easily attached to the light conductor

10447	1 pair of Laryngeal mirrors, 2pcs.



Magnifying glass

- For **ri-scope®** L ophthalmoscope
- Separate magnifying glass with 5x magnification

10965	Separate magnifying glass for **ri-scope®** L ophthalmoscope



Swivel magnifying glass

- For **ri-scope®** F.O. nasal speculum
- 2.5x magnification
- With pin connector

11564	Swivel magnifying glass for **ri-scope®** F.O. nasal speculum



Bulb for pneumatic otoscopy

- Compatible with all **Riester** otoscopes
- For diagnosing middle ear disorders
- Enables optimal eardrum examination through air pressure
- No connector needed – pneumatic bulb links directly via integrated port on the side of the instrument head

10960	Bulb for pneumatic test

906170 Rev. B 2026-08-03

Please visit **riester.de** for more information

Rudolf Riester GmbH P.O. Box 35 | Bruckstraße 31 | 72417 Jungingen | Germany
Tel.: (+49) 7477-9270-0 | Fax.: (+49) 7477-9270-70 | info@riester.de | www.riester.de

Riester