AR-EL TIBBI ÜRÜNLER

2061-2 ENT CHAIR (SWIVEL)

with 2 Motors Electrically adjustable height and backrest Mechanical adjustable swivel





Techinical Details:

Seat Height **:** 590mm - 990 mm
Chair Lenght **:** 1900mm, (without headrest: 1600mm)
Chair Width **:** 650mm, (without armrests: 550mm)
Chair Base **:** 500mm x 700mm
Seat Section **:** 550mm x 450 mm
Back Section **:** 500mm x 540mm
Head Section **:** 210mm x 210mm
Leg Section : 510mm x 260 mm
Foot Step **:** 435mm x 250mm
Motor **:** 2 pieces
Lifting Capacity **:** 200 kg
Command Module **:** Hand swtich

Functional Details:

Rotation **:** 180° left and 180° right and locable.
Head Section **:** Detachable , Adjustable front and back, adjustabe height between 0 -150mm.
Back Section **:** 0° ~ +85°
Foot step **:** Foldable

Packing Details:

Net weight **:** 110 Kg
Gross weight **:** 135 Kg
Packing Dimensions **:** 152x77x99cm

Warranty : 2 years