

VACUUM MATTRESS




34089

• 34089 MAT PLUS VACUUM MATTRESS
PVC, X-ray translucent mattress provided with 4 immobilization straps and 12 transport holes.
Suitable for any kind of emergency situation.
Size: 210x90 cm. Weight: 11.2 kg. Loading capacity: 150 kg.
Supplied with nylon bag size 86x43x17 cm, pump not included.



34088

• 34088 MAT VACUUM MATTRESS
PVC, X-ray translucent mattress with 6 transport handles, perfectly fits the patient body. Sturdy material and careful manufacture make it suitable for Civil Defence, Fire Brigades, Industry and First-aid. Size: 200x75 cm. Weight: 8.8 kg.
Loading capacity: 150 kg.
Supplied with nylon bag size 86x43x15 cm, pump not included.



34091

• 34091 MAT ORANGE VACUUM MATTRESS
X-ray translucent, comfortable vacuum mattress with 8 transport handles and 4 immobilization straps.
Size: 210x90 cm. Weight: 9 kg. Loading capacity: 150 kg.
Supplied with nylon bag size 86x43x15 cm, pump not included.

• 34095 ALUMINIUM PUMP - 0.75 l for vacuum mattress
Rotating handle. Supplied with hose and fitting. See page 395.



Carrying bag included

SEMIRIGID SPLINT KIT



34119

• 34119 SEMIRIGID SPLINT KIT
Kit of 3 splints (arm, wrist, leg) in nylon carrying bag.
All splints are made of polyester + PVC, with lightweight stays providing high strength support.
Adjustable velcro for easy use.
Multilingual instructions: GB, FR, IT, RO.

EXTRICATION KITS


• 34027 SPINAL IMMOBILIZER TYPE A
• 34028 SPINAL IMMOBILIZER TYPE B
GIMA presents 2 different shape of spinal immobilizers, both are ideal for the immobilization and setting free of patients from emergency situations.
They allow to fix vertically and horizontally the head and spine of the patient to limit any injuries caused by wrong transport.
X-ray translucent.
Supplied with different coloured buckles and belt for easy use, padded cushion velcro straps and waterproof bag.



34028 34027

Provide horizontal flexibility and vertical rigidity support

BODY BAGS





• 34116 BODY BAG - PVC - black
Made of PVC (Inner transparent PVC, over black), provided with two long slide fasteners.
Size: 210x90 cm
Loading capacity: 150 kg.

• 34117 BODY BAG - vinyl - blue
Made of vinyl covered nylon conform to standard EN 1865. Provided with a long slide fastener.
Size: 210x90 cm.
Loading capacity: 150 kg.




• 34124 BODY BAG - PVC - black
Made of PVC with one "U" zipper and 4 handles.
Size: 210x90 cm;
0.02 cm thickness
Loading capacity: 150 kg

• 34118 BODY BAG - vinyl - white
Made of vinyl with one "U" zipper and 4 handles.
Size: 210x90 cm;
0.02 cm thickness
Loading capacity: 150 kg



CARRYING COT

34090



With 8 handles for a secure patient transport

• 34090 CARRYING COT
Made of polyester, PVC coated.
Anti-putridity, bacterial and fungus proof, waterproof, easy to wash and disinfect, X-ray translucent material.
Size 190x75 cm.
Loading capacity: 150 kg.
Certified EN1865-1-2010 + A1 2015.
Supplied in a carrying bag with multilingual manual: GB, FR, IT, ES, PT, DE