PATIENT ANTI-HANDLING SUIT

The suit prevents the patient from handling himself during night or day. Back closure with double slider YKK zipper and final reinforcement. The zip is internally protected by fabric, which ensures no contact with patient's skin. The zip also allows the operator to intervene easily on the patient, for hygiene. Long sleeve and leg with elastic cuffs and ankles; double reinforced crotch. Designed for elderly people with incontinence problems associated with neurodegenerative diseases, clothing intolerance or mental disabilities. Made of jersey: 97% cotton and 3% lycra. Washable up to 60°C. Weight: 160 gr/m².



PATIENT ANTI-HANDLING SUIT		
GIMA code	Unisex - green	Size Height x chest
36736	Anti-handling suit - S	150x120 cm
36737	Anti-handling suit - M	152x124 cm
36738	Anti-handling suit - L	154x128 cm
36739	Anti-handling suit - XL	156x134 cm

TOPICAL CREAMS - BODY & FACE

• 28349 TOPICAL CREAM 150 ml - ARNICA AND DEVIL'S CLAW - body
• 28350 TOPICAL CREAM 150 ml - HYALURONIC ACID - face and body
Creams for specific local applications. For external use only.
Arnica and Devil's claw cream: provide relief and comfort. It alleviates discomfort caused by sports and improves skin elasticity in stressed areas. Ideal for anti-fatigue and pre-sport warm-up massages. With its protective and refreshing action, it promotes skin microcirculation and keeps skin toned, radiant and hydrated.
Hyaluronic Acid cream: soothing and calming effect on reddened skin or after sun exposure. It protects the skin from external agents and combats aging. It maintains and improves elasticity and tone, giving firmness and density to sagging and dull skin. It rebalances the hydrolipidic film, keeping skin deeply hydrated and evening out its colour.
Label: GB, IT.



MADE IN ITALY

28349 28350

WIDE SCREEN THERMOMETER

10 SECONDS GREEN/ORANGE/RED BACKLIGHT

Green: lower than 37.4 °C
Orange: between 37.5 and 38.5 °C
Red: between 38.6 °C and 42.9 °C

Multilingual manual and box: GB, FR, IT, ES, PT, SE, HU, RO




Battery 3 V button size (CR2032)

• 25555 3-COLOUR WIDE SCREEN DIGITAL THERMOMETER °C/°F (switchable) - blister
Flexible tip, water-proof, red/orange/green backlight after measurement, wide display 27x20 mm. Low battery indicator

7 DAYS PILL BOX

• 41510 - 7 DAYS PILL BOX blue translucent
With 7 detachable 4 compartments small boxes (morning, noon, evening and night). Translucent top enables the user to see if medication has been taken.
Size: 9.6x2x1.7 cm



It comes with one set of labels for individual marking of compartments in 24. languages: GB, FR, IT, ES, PT, DE, NL,PL, DK, SE, FI, NO, CZ, SK, HU, HR, SI, RO, BG, GR, LT, LV, EE, Arabic.

GENTLE TEMP 730

• 25612 OMRON GENTLE TEMP 730 NON-CONTACT THERMOMETER - MC-F300-E
Gentle Temp 730 is a no touch infrared thermometer using advanced distance sensor technology for fast, accurate and reliable readings of forehead temperature.
It delivers fast and contactless measurement in just one second.
It automatically saves the last 30 readings, making it easy to track temperature trends.
The device can also be switched to surface measurement mode, which shows the actual and unadjusted temperature of an object or liquids. Reading scale °C/°F is switchable.
Buzzer-off and backlight functions.
Ideal for newborns, children and adults.
Supplied with 2xAAA alkaline batteries, sensor cover and multilingual manual: GB, FR, IT, ES, DE, NL, PL, TR, Arabic.



25612

TECHNICAL SPECIFICATIONS	
Sensor:	thermopile
Temperature display:	3-digit, °C display in 0.1 degree increments 4-digit, °F display in 0.1 degree increments
Measuring range:	forehead: +34.0 to +43.0 °C (+93.2 to +109.4 °F), ±0.2 °C (±0.4 °F) surface: 0 to +90.0 °C (+32.0 to +194.0 °F)
Measuring time:	1 second
Memory:	30 readings
Power supply:	DC 3.0 V (2xAAA alkaline batteries)
Battery life:	approx. 1100 measurements
Protection:	type BF, internally powered ME equipment, IP22
Size - weight:	157x27.5xh 37 mm - 102 g

BREAST PUMP - WEARABLE

• 41710 WEARABLE ELECTRIC BREAST PUMP
Designed to offer moms comfort, efficiency, and safety. Thanks to its hands-free, compact design, it fits perfectly inside most nursing bras.
- 2 modes and 9 adjustable vacuum levels
- huge capacity of 180 ml
- removable compartments for sterilization
- quick and gentle suction
- made of food-grade silicone

Supplied with breast shield, silicone diaphragm, linker, valve, milk container, USB charging cable, bra adjustment buckle, multilingual manual (GB, FR, IT, ES).

1.56" LED screen featuring a built-in time display



Transparent design, easily check milk flow and volume

TECHNICAL SPECIFICATIONS	
Breast shield:	24 mm
Capacity:	180 ml
Screen size:	35x18.6 mm
Power supply:	DC 3.7 V - rechargeable Li-battery 1000 mAh
Battery life:	up to 1.8 hours
Size - weight:	110x73.8x130.3 mm - 225 g