VISOR HEAD MASKS

VISOR HEAD MASKS
Great features
- faceshield with 90° rotation
- inner and outer surface treated against mist
- superior coverage to protect against droplets and splashes

Superior materials
- elastic band in anallergic material
- rotation systems in anodized aluminium
- soft support in anallergic interchangeable foam
- faceshield in heat-formed acetate with anti-fog treatment

Easy and comfortable use
- fitting headbands adjustable by velcro fasteners
- usable over ophthalmic glasses

GIMA code	VISOR HEAD	Minimum order
25656	VISOR HEAD MASK	1
25668	VISOR SHIELD	1
25669	SPARE SHIELD for 25668	box of 5




TETI PROTECTIVE SHIELDS

- **30863 TETI PROTECTIVE SHIELD - white - 1 shield**
- **30860 TETI PROTECTIVE SHIELD - white - 2 shields**

Comfortable, ultra light stylish protective shields. They provide high protection from squirts and sprays.
Folding shield behind the nape.
Adjustable size, fits majority of healthcare professionals.
Extra wide visor allows to fit over glasses and most of magnification loupes.
Supplied with 2 face transparent shields and 2 padding with velcro (30860) or 1 shield + 1 pack (30863)
Certificated EN166. Made in Italy.

- **30875 SET OF 5 TRANSPARENT FACE SHIELDS - spare**



OP-DOP VISOR SHIELD "PROTECTOR" WITH 3+1 SHIELDS

- **25642 VISOR SHIELD "PROTECTOR" - white**

The ultra-light face shield that combines function with style:
protection: a waterproof plastic barrier will protect you
comfort: fast and easy to wear and remove
convenience: reusable and replaceable shield
optical quality: superior optical quality.
Medium adjustable size: fits over 90% of healthcare professionals. Fits over glasses.
Shielded air vent eliminates fogging.
Supplied with 3 face transparent shields + 1 head shield.

- **25643 REPLACEMENT SHIELD KIT**


Includes 2 face transparent shields + 1 head shield + 1 foam bumper





Quality product: made in USA

GOGGLES KIT



- **25646 GOGGLES KIT - fog resistant**

Practical kit composed of 5 supports in different colours and 20 antifog shields. Usable over ophthalmic glasses.

- **25647 SPARE SHIELDS - box of 20**

HIGH PROTECTION MASK 6x1

Instructions in: GB, FR, IT, DE, ES, PT, NL, DK, NO, SE, FI, GR, PL, CZ, SI, HR, HU, RO, BG, RU, TR, SA, JP.



- **25265 HIGH PROTECTION MASK 6x1 - fog resistant, anti-scratch**
- wide lens for a field of view without distortion
- temples adjustable in length and inclination with SoftPad technology
- sophisticated indirect ventilation system
- soft internal rubber for comfort without local pressure
- certified according to EN166 point 3 (droplets and splashes) and point 4 (large dust particles)
- supplied with neck cord